DIAQUICK DOA Cassettes

for human urine samples

	REF	Content
DIAQUICK AMP Cassette	Z99004CE	- 30 Tests (30x REF Z99004B)
DIAQUICK BAR Cassette	Z99006CE	- 30 Tests (30x REF Z99006B)
DIAQUICK BUP Cassette	Z04560CE	- 30 Tests (30x REF Z04560B)
DIAQUICK BZO Cassette	Z99001CE	- 30 Tests (30x REF Z99001B)
DIAQUICK COC Cassette	Z99003CE	- 30 Tests (30x REF Z99003B)
DIAQUICK ETG Cassette	Z15102CE	- 30 Tests (30x REF Z15102B)
DIAQUICK FYL Cassette	Z09640CE	- 10 Tests (10x REF Z09640B)
DIAQUICK KET Cassette	Z09641CE	- 10 Tests (10x REF Z09641B)
DIAQUICK MDMA Cassette	Z04570CE	- 30 Tests (30x REF Z04570B)
DIAQUICK MET Cassette	Z99500CE	- 30 Tests (30x REF Z99500B)
DIAQUICK MOP Cassette	Z99005CE	- 30 Tests (30x REF Z99005B)
DIAQUICK MTD Cassette	Z99550CE	- 30 Tests (30x REF Z99550B)
DIAQUICK OPI Cassette	Z05011CE	- 30 Tests (30x REF Z05011B)
DIAQUICK TCA Cassette	Z03040CE	- 30 Tests (30x REF Z03040B)
DIAQUICK THC Cassette	Z99002CE	- 30 Tests (30x REF Z99002B)
DIAQUICK TRA Cassette	Z10414CE	- 30 Tests (30x REF Z10414B)
DIAQUICK Spice Cassette	Z13630CE	- 30 Tests (30x REF Z13630B)

All tests are individually packed and contain a disposable plastic pipette.
 All products contain a package insert.

For in vitro diagnostic use only. For use by medical professionals only.
 For diagnosis and therapeutic monitoring only.

INTENDED USE

The DIAQUICK DOA Cassettes (urine) are rapid, lateral flow chromatographic immunoassays for the qualitative detection of the following drugs and their metabolites:

Parameter	Code	Calibrator Substance	Cut-off
Amphetamine	AMP	d-Amphetamine	1 000 ng/mL
Barbiturates	BAR	Secobarbital	300 ng/mL
Buprenorphine	BUP	Buprenorphine	10 ng/mL
Benzodiazepines	BZO	Oxazepam	300 ng/mL
Cocaine	COC	Benzoylcegonine	300 ng/mL
Ethylglucuronide	ETG	Ethyl-β-D-Glucuronide	500 ng/mL
Fentanyl	FYL	Norfentanyl	20 ng/mL
Ketamine	KET	Ketamine	1 000 ng/mL
Ecstasy	MDMA	(±)3,4-Methylenedioxymethamphetamine HCl	500 ng/mL
Methamphetamine	MET	d-Methamphetamine	1 000 ng/mL
Opiate, Morphine, Heroin	MOP	Morphine	300 ng/mL
Methodone	MTD	Methodone	300 ng/mL
Opiate, Morphine, Heroin	OPI	Morphine	2 000 ng/mL
Tricyclic Antidepressants	TCA	Nortriptyline	1 000 ng/mL
Marihuana/Cannabis	THC	11-nor-Δ9-THC-9-COOH	50 ng/mL
Tramadol	TRA	cis-Tramadol	100 ng/mL
Synthetic Marihuana	K2	JWH-018 5-Pentanoic acid	50 ng/mL

These tests will detect other related compounds; please refer to the Analytical Specificity table in this insert. This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are obtained.

TEST PRINCIPLE

The DIAQUICK DOA Cassettes (urine) are immunoassays based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against their respective drug conjugate for binding sites on their specific antibody. During testing, a urine specimen migrates upward by capillary action. A drug, if present in the urine specimen below its cut-off concentration, will not saturate the binding sites of the specific antibody coated on the particles. The antibody coated particles will then be captured by the immobilized drug conjugate and a visible coloured line will show up in the test line region. The coloured line will not form in the test line region if the drug level is above its cut-off concentration because it will saturate all the binding sites of the antibody coated on the particles. A drug-positive urine specimen will not generate a coloured line in the specific test line region of the strip because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a coloured line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

REAGENTS

The test line contains mouse monoclonal antibody-coupled particles and corresponding drug-protein conjugates. A goat antibody is employed in the control line.

WARNINGS AND PRECAUTIONS

- For medical and other in vitro diagnostic use only. Do not use after the expiration date.
- The test cassettes should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test cassettes should be discarded according to federal, state and local regulations.

STORAGE

The DIAQUICK DOA Cassettes can be stored refrigerated or at room temperature (2 – 30 °C). The tests are stable through the expiration date printed on the sealed pouch. The test cassettes must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

SAMPLE COLLECTION AND PREPARATION

The urine must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitations should be centrifuged, filtered or allowed to settle to obtain a clear specimen for testing. Urine specimens may be stored at 2 – 8 °C for up to 48 h prior to testing. For prolonged storage, specimens may be frozen and stored below –20 °C. Frozen specimens should be thawed and mixed before testing.


MATERIALS REQUIRED BUT NOT PROVIDED

- Specimen collection container
- Timer

ASSAY PROCEDURE

Allow the test cassette, urine specimen, and/or controls to equilibrate to room temperature (15 – 30 °C) prior to testing.

- Bring the pouch to room temperature before opening it. Remove the test cassette from the sealed pouch and use it as soon as possible.
- Place the test cassette on a clean and level surface. Hold the pipette vertically and transfer 3 full drops of urine (approx. 120 µL) to the specimen well (S) of the test cassette and start the timer. Avoid trapping air bubbles in the specimen well (S). See the illustration below.
- Wait for the coloured line(s) to appear. Read results at 5 minutes. Do not interpret the result after 10 minutes.

INTERPRETATION OF RESULTS

NEGATIVE: Two lines appear. One coloured line should be in the control line region (C), and another apparent coloured line should be in the test line region (T). This negative result indicates that the drug concentration is below the detectable level.

*NOTE: The shade of colour in the test region (T) may vary, but it should be considered negative whenever there is even a faint coloured line.

POSITIVE: One coloured line appears in the control line region (C). No line appears in the test line region (T). This positive result indicates that the drug concentration exceeds the detectable level.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test cassette. If the problem persists, discontinue using the lot immediately and contact your local distributor.

QUALITY CONTROL

A procedural control is included in the test. A coloured line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

- The DIAQUICK DOA Cassettes provide only a preliminary analytical result. A more specific chemical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.^{1,2}
- It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
- Adulterants, such as bleaching agents in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
- A positive result indicates presence of the drug or its metabolites but does not indicate the level of intoxication, administration route or concentration in urine.
- A negative result may not necessarily indicate drug-free urine. Negative results can be obtained if a drug is present but below the cut-off level of the test.
- The DIAQUICK DOA Cassettes do not distinguish between drugs of abuse and certain medications.
- A positive result might be obtained from certain foods or food supplements.

PERFORMANCE CHARACTERISTICS

ACCURACY

A side-by-side comparison of the DIAQUICK DOA Cassettes and a commercially available rapid drug test was conducted. Testing was performed on approx. 100 specimens previously collected from subjects present for drug screen testing. The agreement was > 99.9 % for all tests.

A side-by-side comparison of the DIAQUICK DOA Cassettes and GC/MS at the cut-off level of the tests was conducted. Testing was performed on 250 specimens previously collected from subjects present for drug screen testing. The following results were tabulated:

	% Agreement with GC/MS		
	Positive Agreement	Negative Agreement	Total Results
AMP	98,1 %	97,9 %	98,0 %
BAR	96,1 %	98,6 %	97,6 %
BUP	99,1 %	> 99,9 %	99,6 %
BZO	98,4 %	99,2 %	98,8 %
COC	98,2 %	97,8 %	98,0 %
ETG	97,6 %	99,4 %	98,8 %
FYL	98,8 %	99,4 %	99,2 %
KET	97,5 %	98,2 %	98,0 %
MDMA	98,1 %	99,3 %	98,8 %
MET	96,2 %	97,1 %	96,8 %
MOP	95,0 %	95,3 %	95,2 %
MTD	98,9 %	98,8 %	98,8 %
OPI	96,7 %	93,8 %	95,2 %
TCA	94,8 %	91,6 %	92,8 %
THC	97,9 %	98,1 %	98,0 %
TRA	88,2 %	92,4 %	90,8 %
K2	97,5 %	98,2 %	98,0 %

ANALYTICAL SPECIFICITY

The following tables lists the concentration of compounds (ng/mL) that are detected positive

in urine by the DIAQUICK DOA Cassettes at 5 minutes.

AMPHETAMINE	AMP	BARBITURATES	BAR
D,L-Amphetamine sulfate	300	Amobarbital	5 000
L-Amphetamine	25 000	5,5-Diphenylhydantoin	8 000
(±) 3,4-Methylenedioxyamphetamine	500	Allobarbitol	600
Phentermine	800	Barbital	8 000
Maprotiline	50 000	Talbutal	200
Methoxyphenamine	6 000	Butalbutal	8 000
D-Amphetamine	1 000	Phenobarbital	300
BUPRENORPHINE	BUP	Cyclopentobarbital	30 000
Buprenorphine	10	Pentobarbital	8 000
Norbuprenorphine	50	Alphenol	600
Buprenorphine 3-D-Glucuronide	50	Aprobarbital	500
Norbuprenorphine 3-D-Glucuronide	100	Butabarbitol	200
BENZODIAZEPINES	BZO	Butethal	500
Alprazolam	100	Secobarbital	300
a-hydroxyalprazolam	1 500	COCAINE	COC
Bromazepam	900	Benzoylcegonine	300
Chlordiazepoxide	900	Cocaine HCl	200
Clobazam	200	Cocaehtylene	20 000
Clonazepam	500	Ecgonine HCl	30 000
Clorazepate dipotassium	500	ETHYLGLUCURONIDE	ETG
Delorazepam	900	Ethyl-β-D-Glucuronide	500
Desalkylflurazepam	200	Propyl-β-D-Glucuronide	50 000
Diazepam	300	Morphine-3-β-Glucuronide	100 000
Estazolam	6 000	Morphine-6-β-Glucuronide	100 000
Flunitrazepam	200	Glucuronic Acid	100 000
(±) Lorazepam	3 000	Ethanol	100 000
RS-Lorazepam glucuronide	200	Methanol	100 000
Midazolam	6 000	FENTANYL	FYL
Nitrazepam	200	Alfentanil	600 000
Norchlordiazepoxide	100	Fenfluramine	50 000
Nordiazepam	900	Norfentanyl	20
Oxazepam	300	Busporine	15 000
Temazepam	100	Fentanyl	100
Triazolam	3 000	Sufentanyl	50 000
KETAMINE	KET	ECSTASY	MDMA
Ketamine	1 000	(±) 3,4-Methylenedioxyamphetamine HCl	500
Benzphetamine	25 000	(±) 3,4-Methylenedioxyamphetamine HCl (MDA)	3 000
(+) Chlorpheniramine	25 000	3,4-Methylenedioxyethyl-amphetamine (MDE)	300
Clonidine	100 000	METHAMPHETAMINE	MET
Dextromethorphan	2 000	µ-Hydroxymethamphetamine	25 000
Disopyramide	25 000	D-Methamphetamine	1 000
EDDP	50 000	L-Methamphetamine	20 000
Mephentermine	25 000	(±)-3,4-Methylenedioxyamphetamine	12 500
(1R, 2S) - (-)-Ephedrine	100 000	Mephentermine	50 000
4-Hydroxyphencyclidine	50 000	MORPHINE	MOP
Levorphanol	50 000	Codeine	200
MDE	50 000	Ethylmorphine	6 000
Tetrahydrozoline	500	Hydrocodone	50 000
d-Methamphetamine	50 000	Hydromorphone	3 000
l-Methamphetamine	50 000	Levorphanol	1 500
Methoxyphenamine	25 000	6-Monoacetylmorphine	300
(+)-3,4-Methylenedioxyamphetamine	100 000	Morphine 3-β-D-glucuronide	800
d-Norpropoxyphene	25 000	Morphine	300
Pentazocine	25 000	Norcodeine	6 000
Phencyclidine	25 000	Normorphone	50 000
Promazine	25 000	Oxycodone	30 000
Promethazine	25 000	Oxymorphone	50 000
Thioridazine	50 000	Procaine	15 000
Meperidine	25 000	Thebaine	6 000
CANNABIS	THC	TRICYCLIC ANTIDEPRESSANTS	TCA
Cannabinol	35 000	Nortriptyline	1 000
11-nor-Δ ⁸ -THC-9 COOH	30	Nordoxepine	500
11-nor-Δ⁹-THC-9 COOH	50	Trimipramine	3 000
Δ ⁸ -THC	17 000	Amitriptyline	1 500
Δ ⁹ -THC	17 000	Promazine	3 000
SPICE	K2	Desipramine	200
JWH-018 5-Pentanoic acid metabolite	50	Cyclobenzaprine	2 000
JWH-073 4-butanoic acid metabolite	50	Imipramine	400
JWH-018 4-Hydroxypentyl metabolite	400	Clomipramine	50 000
JWH-018 5-Hydroxypentyl metabolite	500	Doxepine	2 000
JWH-073 4-Hydroxybutyl metabolite	500	Maprotiline	2 000
JWH-073 N-(3-hydroxypentyl) metabolite	8 000	Promethazine	50 000
JWH-018 N-(4-hydroxypentyl) metabolite	10 000	Perphenazine	50 000
MAM2201 N-Pentanoic metabolite	300	Dithiaden	10 000
JWH-122 N-(4-hydroxypentyl) metabolite	2 000	METHADONE	MTD
JWH-018 N-Pentanoic metabolite	150	Methadone	300
JWH-073 N-(2-hydroxybutyl) metabolite	5 000	Doxylamine	100 000
JWH-018 N-(5-hydroxypentyl) metabolite	5 000	Cis-tramadol	300 000
JWH-019 5-hydroxypentyl metabolite	10 000	OPIATES	OPI
JWH-019	10 000	Codeine	2 000
JWH-122 N-(5-hydroxypentyl) metabolite	5 000	Ethylmorphine	3 000
JWH-398 N-Pentanoic acid metabolite	500	Hydrocodone	50 000
JWH-200 6-hydroxyindole metabolite	15 000	Hydromorphone	15 000
JWH-210 N-Pentanoic acid metabolite	1 000	Levorphanol	25 000
RCS4 N-5-Carboxypentyl metabolite	1 000	6-Monoacetylmorphine	3 000
JWH-073 4-Pentanoic acid metabolite	10 000	Morphine 3-β-D-glucuronide	2 000
TRAMADOL	TRA	Morphine	2 000
n-Desmethyl-cis-tramadol	200	Norcodeine	25 000
Cis-tramadol	100	Normorphone	50 000
Procyclidine	100 000	Oxycodone	25 000
o-Desmethyl-cis-tramadol	10 000	Oxymorphone	25 000
Phencyclidine	100 000	Procaine	50 000
d,l-O-Desmethyl venlafaxine	50 000	Thebaine	25 000

CROSS-REACTIVITY


A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free or drug positive urine. The following compounds did not show a cross-reactivity when tested with the DIAQUICK DOA Cassettes at a concentration of 100 µg/mL.

Non Cross-Reacting Compounds:

Acetophenetidin	Cortisone	Zomepirac	d-Pseudoephedrine
N-Acetylprocainamide	Creatinine	Ketoprofen	Quinidine
Acetylsalicylic acid	Deoxycorticosterone	Labeltalol	Quinine
Aminopyrine	Dextromethorphan	Loperamide	Salicylic acid
Amoxicillin	Diclofenac	Meprobamate	Serotonin
Ampicillin	Diffunisal	Methoxyphenamine	Sulfamethazine
l-Ascorbic acid	Digoxin	Methylphenidate	Sulindac
Apomorphine	Diphenhydramine	Nalidixic acid	Tetracycline
Aspartame	Ethyl-p-aminobenzoate	Naproxen	Tetrahydrocortisone,
Atropine	β-Estradiol	Niacinamide	3-acetate
Benzilic acid	Estrone-3-sulfate	Nifedipine	Tetrahydrocortisone
Benzoic acid	Erythromycin	Norethindrone	Tetrahydrozoline
Bilirubin	Fenoprofen	Noscapine	Thiamine
d,l-Brompheniramine	Furosemide	d,l-Octopamine	Thioridazine
Caffeine	Gentisic acid	Oxalic acid	d,l-Tyrosine
Cannabidiol	Hemoglobin	Oxolinic acid	Tolbutamide
Chloral hydrate	Hydralazine	Oxymetazoline	Triamterene
Chloramphenicol	Hydrochlorothiazide	Papaverine	Trifluoperazine
Chlorothiazide	Hydrocortisone	Penicillin-G	Trimethoprim
d,l-Chlorpheniramine	o-Hydroxyhippuric acid	Perphenazine	d,l-Tryptophan
Chlorpromazine	3-Hydroxytyramine	Phenelzine	Uric acid
Cholesterol	d,l-Isoproterenol	Prednisone	Verapamil
Clonidine	Isoxsuprine	d,l-Propranolol	

REFERENCES

- Baselt, RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488
- Tietz NW. Textbook of Clinical Chemistry. W.B. Saunders Company. 1986; 1735
- Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986.

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