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Instruction for use

A solid-phase enzyme immunoassay for the quantitative determination of IgG antibodies to deamidated gliadin peptides (DGP) in human blood serum or plasma

1. INTENDED USE

A solid-phase enzyme immunoassay for the quantitative determination of IgG antibodies to deamidated gliadin peptides (DGP) in blood serum or plasma.

This kit is designed for measurement of IgG antibodies to deamidated gliadin peptides (DGP) in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 41 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Celiac disease (CD) is a complex autoimmune pathology characterized by intolerance to a major reserve protein of rye and wheat grains – gluten. Prevalence of this disease varies greatly in different geographic regions and may come to 1:300. If a patient with CD consumes food products containing gluten, a crypt hyperplasia accompanied by total or partial atrophy of intestinal villi develops; at the same time, the risk of development of oncological pathology of the intestine rises significantly.

Celiac diagnostic procedures include determination of anti-tissue transglutaminase auto antibodies (IgA and IgG), anti-gliadin auto antibodies (IgA and IgG), anti-endomysium antibodies (EMA) and biopsy of small intestine as a main confirmation test. Auto antibodies against deamidated gliadin peptides (DGP) seem to be more specific marker in addition to determination of anti-gliadin auto antibodies. During pathogenesis of celiac disease gliadin undergoes deamidation by tissue transglutaminase in the lamina propria. Deamidated gliadin peptides bind to HLA molecules on the surface of antigen presenting cell and are recognized by T-helper cells. It causes anti-DGP antibodies appearance. Therefore determination of anti-DGP auto antibodies may support diagnostic of celiac disease especially as addition marker to anti – total gliadin auto antibodies determination.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by the antigen. Antibodies from the specimen bind coated antigen on the microwell surface. Unbound material is removed by washing procedure. Second antibodies directed towards species specific Ig, labelled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

4.1. For professional use only.

4.2. This kit is intended for in vitro diagnostic use only.

4.3. INFECTION HAZARD: There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

4.4. Avoid contact with stop solution containing 5,0 % H_2SO_4 . It may cause skin irritation and burns.

4.5. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

4.6. Do not use the kit beyond the expiration date.

4.7. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

4.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

4.9. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

4.10. Do not mix reagents from different lots.

4.11. Replace caps on reagents immediately. Do not swap caps.

4.12. Do not pipette reagents by mouth.

4.13. Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

4.14. Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

4.15. The Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	Anti-DGP IgG EIA strips, 8x12 wells	1	pcs		until exp.date
2	Calibrator set, 1.5 ml each. The set contains 6 calibrators: 0; 5; 10; 25; 50; 100 IU/ml	6	pcs	red (C1 - colourless)	2 months
3	Control serum (1.5 ml)	1	pcs	colourless	2 months
4	Conjugate, 11 ml	1	pcs	red	until exp.date
5	EIA buffer 50 ml	1	pcs	blue	until exp.date
6	Substrate solution, 11 ml	1	pcs	colourless	until exp.date
7	Washing solution concentrate 21X; 22 ml	1	pcs	colourless	Concentrate - until exp.date Diluted washing solution - 1 month at 2-8 °C or 5 days at RT
8	Stop solution, 1.1 ml	1	pcs	colourless	until exp.date
9	Plate sealing tape	2	pcs		N/A
10	Instruction Anti-DGP IgG EIA	1	pcs		N/A
11	QC data sheet Anti-DGP IgG EIA	1	pcs		N/A

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 5–250 µl;
- Dry thermostat for +37 °C ±0.1 °C
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0.

5.3. Storage and stability of the Kit

Store the whole kit at +2...+8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells **TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED)** to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at +2...+8 °C before testing. Freezing/thawing should be avoided.

7. TEST PROCEDURE**7.1.** Reagent Preparation

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 21X by 21 dilutions in distilled water.

7.2. Procedural Note:

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Alternative units:

7.5. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 14 wells for the calibrators CAL 1-6 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Dilute all samples using buffer DIL SPE (EIA buffer) 101 fold. Do not dilute control sample and calibrators.
3	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using DIL SPE (EIA buffer). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
4	Pipet 100 µl of calibrators CAL 1-6 and unknown samples into the wells.
5	Incubate 30 minutes at +37 °C .
6	Prepare washing solution by 21x dilution of washing solution concentrate BUF WASH 21X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times
7	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
8	Incubate 30 minutes at +37 °C .
9	Wash the strips 5 times.
10	Dispense 100 µl of SUBS TMB into the wells
11	Incubate 10-20 minutes at +18...+25 °C
12	Dispense 100 µl of STOP into the wells.
13	Measure OD (optical density) at 450 nm.
14	Set photometer blank on first calibrator
15	Apply point-by-point method for data reduction.

7.6. Sample processing

Material type	Notes on material collection, storage and handling	Sample dilution example	Sample into the well, µl	Calculation factor
blood serum or plasma	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	5 µl of sample + 500 µl of diluent	1.00	1

8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of *GLP* (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

9. CALCULATION OF RESULTS


9.1. Calculate the mean absorbance values (OD450) for each pair of calibrators and samples.

9.2. Plot a calibration curve on graph paper: OD versus IgG antibodies to deamidated gliadin peptides (DGP) concentration.

9.3. Determine the corresponding concentration of IgG antibodies to deamidated gliadin peptides (DGP) in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

9.4. Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

Calibrators	Value	Absorbance Units (450 nm)
CAL 1	0 IU/ml	0.07
CAL 2	5 IU/ml	0.18
CAL 3	10 IU/ml	0.54
CAL 4	25 IU/ml	0.94
CAL 5	50 IU/ml	1.67
CAL 6	100 IU/ml	2.32

10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for Anti-DGP IgG.

Sex, age	Units, U/ml	
	Lower limit	Upper limit
Healthy donors	-	10

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical sensitivity

Sensitivity of the assay was assessed as being 5 IU/ml.

11.2. Linearity

Linearity was checked by assaying dilution series of 5 samples with different IgG antibodies to deamidated gliadin peptides (DGP) concentrations. Linearity percentages obtained ranged within 90 to 110 %.

11.3. Recovery

Recovery was estimated by assaying 5 mixed samples with known IgG antibodies to deamidated gliadin peptides (DGP) concentrations. The recovery percentages ranged from 90 to 110 %.