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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF TSH IN HUMAN BLOOD SERUM OR PLASMA

1. INTENDED USE

A solid-phase enzyme immunoassay for the quantitative determination of TSH in blood serum or plasma.

This kit is designed for measurement of TSH in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 39 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Thyroid stimulating hormone (TSH) is a glycoprotein with molecular weight ca.30 kDa which is secreted by hypophysis. A molecule of TSH consists of two noncovalently bound subunits: α - and β -HCG. β -subunit determines biological activity and immunological specificity of TSH.

TSH stimulates thyroid gland to secrete thyroid hormones. TSH secretion in hypophysis is controlled by a negative feedback regulation by thyroid hormones. TSH secretion is subject to circadian rhythms with highest levels seen early in the morning (6 a.m.). Changes of TSH blood level during a day are not significant; nevertheless, if the results do not correspond with clinical status and other laboratory data, it is recommended to take and test another blood sample.

Determination of TSH level in serum is recommended in the following states and conditions:

1. Diagnostics of dysfunction of the thyroid gland;
2. Hypothyroidism (TSH level is increased. The diagnosis is confirmed by low concentrations of total and free T4 and T3. In mild subclinical forms when T4 and T3 levels are within normal ranges, determination of TSH concentration is critical);
3. Hyperthyroidism (synthesis and secretion of TSH are inhibited); monitoring of replacement therapy;
4. Screening for inherited hypothyroidism (on day 5 after birth TSH level in blood is determined). TSH level is elevated just after birth but it comes within the normal range in several days (both in boys and in girls).

Serum TSH level is elevated during pregnancy, after physical stress, in individuals with lowered blood pressure and lowered temperature. Secretion of TSH is inhibited by Cortisol and Growth hormone. Low TSH levels are often seen in elderly people, in patients with chronic renal insufficiency, liver cirrhosis, in retardation of sexual development, in secondary amenorrhea, Cushing syndrome, acromegaly.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to β chain human TSH-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Unbound material is removed by washing procedure. Second antibodies - murine monoclonal to (Fab2)-fragment of β chain human TSH, labelled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

4.1. For professional use only.

4.2. This kit is intended for in vitro diagnostic use only.

4.3. INFECTION HAZARD: There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

4.4. Avoid contact with stop solution containing 5,0% H₂SO₄. It may cause skin irritation and burns.

4.5. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

4.6. Do not use the kit beyond the expiration date.

4.7. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

4.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

4.9. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

4.10. Do not mix reagents from different lots.

4.11. Replace caps on reagents immediately. Do not swap caps.

4.12. Do not pipette reagents by mouth.

4.13. Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

4.14. Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

4.15. The Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	TSH plus EIA strips, 8x12 wells	1	pcs		until exp.date
2	Calibrator set, 1 ml each, zero calibrator C1 - 2 ml The set contains 7 callibrators: 0; 0,05; 0,2; 1; 5; 15, 50 mIU/l	7	pcs	red (C1 - colourless)	2 months
3	CONTROL _L , CONTROL ₂	2	pcs	colourless and colourless	2 months
4	CONJ HRP	1	pcs	bright red	until exp.date
5	DIL	1	pcs	blue	until exp.date
6	SUBS TMB	1	pcs	colourless	until exp.date
7	BUF WASH 21X	1	pcs	colourless	Concentrate - until exp.date Diluted washing solution - 1 month at 2-8 °C or 5 days at RT
8	STOP	1	pcs	colourless	until exp.date
9	I003	2	pcs		N/A
10	K201AI	1	pcs		N/A
11	K201AQ	1	pcs		N/A

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 25–250 µl;
- Dry thermostat for 37 °C ±0,1 °C;
- Calibrated microplate photometer with 450 nm and 492 nm wavelength and OD measuring range 0-3.0.

5.3. Storage and stability of the Kit

Store the whole kit at 2 to 8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells **TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED)** to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at 2-8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

7. TEST PROCEDURE**7.1. Reagent Preparation**

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18 to +25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 21X by 21 dilutions in distilled water.

7.2. Procedural Note:

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Assay procedure

Sensitivity 0.025 mIU/l

1	Put the desired number of microstrips into the frame; allocate 18 wells for the calibrators CAL 1 - 7 and control samples CONTROL, CONTROL 2 and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using (zero calibrator). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
3	Pipet 100 µl of EIA buffer into each well.
4	Pipet 100 µl of calibrators CAL 1 - 7, control samples CONTROL, CONTROL 2 and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
5	Incubate 60 minutes at 37 °C .
6	Prepare washing solution by 21x dilution of washing solution concentrate BUF WASH 21X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times
7	Dispense 200 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
8	Incubate 60 minutes at 37 °C .
9	Wash the strips 5 times.
10	Dispense 200 µl of SUBS TMB into the wells
11	Incubate 10-20 minutes at 18-25 °C
12	Dispense 100 µl of STOP into the wells.
13	Measure OD (optical density) at 450 nm.
14	Set photometer blank on first calibrator
15	Apply point-by-point method for data reduction in the range of 0-5 mIU/l. ATTENTION: 1.5 and 50 mIU/l points should be omitted in data reduction by standard run mode! Repeat OD measurement at 492 nm. Set photometer blank on first calibrator. Apply point-by-point method for data reduction in the range of 0-50 mIU/l. ATTENTION: 0.05 mIU/l point should be omitted in data reduction by standard run mode!

Fast run mode
Sensitivity 0.1 mIU/l

1	Put the desired number of microstrips into the frame; allocate 18 wells for the calibrators CAL 1 - 7 and control samples CONTROL, CONTROL 2 and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using zero calibrator. Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
3	Pipet 100 µl of EIA buffer into each well.
4	Pipet 100 µl of calibrators CAL 1 - 7, control samples CONTROL, CONTROL 2 and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
5	Incubate 30 minutes at 37 °C.
6	Prepare washing solution by 21x dilutions of washing solution concentrate BUF WASH 21X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times
7	Dispense 200 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
8	Incubate 30 minutes at 37 °C.
9	Wash the strips 5 times.
10	Dispense 200 µl of SUBS TMB into the wells
11	Incubate 10-20 minutes at 18-25 °C
12	Dispense 100 µl of STOP into the wells.
13	Measure OD (optical density) at 450 nm.
14	Set photometer blank on first calibrator
15	Apply point-by-point method for data reduction in the range of 0-15 mIU/l. ATTENTION: 50 mIU/l point should be omitted in data reduction by standard run mode! Repeat OD measurement at 492 nm. Set photometer blank on first calibrator. Apply point-by-point method for data reduction in the range of 0-50 mIU/l. ATTENTION: 0.05 mIU/l point should be omitted in data reduction by standard run mode!

8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

9. CALCULATION OF RESULTS


9.1. Calculate the mean absorbance values (OD450) for each pair of calibrators and samples.

9.2. Plot a calibration curve on graph paper: OD versus TSH concentration.

9.3. Determine the corresponding concentration of TSH in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

9.4. Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

Calibrators	Value	Absorbance Units (450 nm) x1000 (fast run OD450)
CAL 1	0 mIU/l	0
CAL 2	0.05 mIU/l	27
CAL 3	0.2 mIU/l	67
CAL 4	1 mIU/l	225
CAL 5	5 mIU/l	890
CAL 6	15 mIU/l	2079

10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for TSH plus. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units	
	mIU/l	
	Lower limit	Upper limit
Healthy donors	0.3	4.0

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical specificity / Cross reactivity

Analyte	Cross-reactivity, % wt/wt
HCG	<0.1
LH	<0.1
FSH	<0.1

11.2. Analytical sensitivity.

Sensitivity of the assay was assessed as being 0.025 mIU/l.

11.3. Linearity.

Linearity was checked by assaying dilution series of 5 samples with different TSH concentrations. Linearity percentages obtained ranged within 90 to 110%.

11.4. Recovery.

Recovery was estimated by assaying 5 mixed samples with known TSH concentrations. The recovery percentages ranged from 90 to 110%.

12. LITERATURE

1. Ekins R. Methods for measurement of free thyroid hormones. In: Free thyroid hormones. Amsterdam: Expecta Medica; 1979; p. 72-92.
2. Tietz, N., Fundamentals of Clinical Chemistry, W.B. Saunders Co., Philadelphia: 791 and 844 (1976).
3. Soos, M., Taylor, S.J., Gard, T., and Siddle, K.A. Rapid Sensitive Two-Site Immunometric Assay for TSH Using Monoclonal Antibodies: Investigation of Factors Affecting Optimisation. J. of Immunological Methods 73, 237-249 (1984).
4. Musto, J.D., Pizzolante, J.M., Chesarone, V.P. A Comment of Thyrotropin Measurement and Evaluation. Clin. Chem. 30, 329-330 (1984). Opinion.
5. Burger, H. G., Patel, Y. C., Thyrotropin releasing hormone -TSH Clinic. Endocrinol. and Metab., 6, 831-00(1977).
6. Ezrin, C., The Thyroid, S. C. Werner and S. H. Ingbar (eds.), Harper and Row, Hagerstown, MD, 9, 174-178 (1978).
7. Engall, E., Methods in Enzymology, Volume 70, Van Vunakis, H. and Langone, J. J. (eds.), Academic Press, New York, 419-492(1980).