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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF PROLACTIN IN HUMAN BLOOD SERUM OR PLASMA

1. INTENDED USE

A solid-phase enzyme immunoassay for the quantitative determination of Prolactin in blood serum or plasma.

This kit is designed for measurement of Prolactin in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 42 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Prolactin is a 198 aminoacids polypeptide with molecular weight ca. 22.5 kDa which is secreted by eosinophil cells of hypophysis.

Hyperplasia and adenomas of hypophysis are the main causes of infertility. Functional changes in regulation of reproductory function are also caused by alterations in secretion of hormones of hypophysis. One of markers of such alterations is changes in Prolactin secretion. That is why WHO recommended to use determination of Prolactin level as a screening test in primary laboratory investigation of couples claiming for infertility.

In women Prolactin level remains low before menarche and elevates during puberty. During this period, Prolactin stimulates development of mammary glands. Prolactin level changes during menstrual cycle with elevations up to 900 mIU/l seen during periovulatory period and the second stage of luteinic phase. That is why it is recommended to evaluate Prolactin level during the first stage of the cycle. Besides, physiological hyperprolactinemia is seen in stress conditions and after physical exercises.

Prolactin secretion is subject to circadian rhythms with maximal levels found during the night (3-7 fold higher than during the day). That is why time of sampling is extremely important.

Elevated Prolactin levels are seen in Prolactin-producing tumours of hypophysis, idiopathic hyperprolactinemias (symptoms: in women – alteration of menstrual cycle, in men – impotence), hypofunction of the thyroid gland, renal insufficiency, after intake of phenothiazine derivatives, haloperidol, estrogens, oral contraceptives, histamine preparations, opiates, in hypoglycemia caused by insuline intake.

Low Prolactin levels are found after surgical resection of hypophysis, after X-ray therapy, after bromocriptin therapy, after intake of T4.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to human Prolactin-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Second antibodies – murine monoclonal to human Prolactin, labelled with peroxidase enzyme, are then added into the microwells. After washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

4.1. For professional use only.

4.2. This kit is intended for in vitro diagnostic use only.

4.3. INFECTION HAZARD: There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

4.4. Avoid contact with stop solution containing 5,0% H₂SO₄. It may cause skin irritation and burns.

4.5. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

4.6. Do not use the kit beyond the expiration date.

4.7. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

4.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

4.9. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

4.10. Do not mix reagents from different lots.

4.11. Replace caps on reagents immediately. Do not swap caps.

4.12. Do not pipette reagents by mouth.

4.13. Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

4.14. Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

4.15. The Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	SORB MTP Prolactin EIA strips, 8x12 wells	1	pcs		until exp. date
2	CAL 1-5 Calibrator set, 0.8 ml each, zero calibrator C1 – 2 ml. The set contains 5 calibrators: 0; 100; 200; 1000; 2000 mIU/l	5	pcs	red (C1 – colourless)	2 months
3	CONTROL Control serum (0.8 ml)	1	pcs	colourless	2 months
4	CONJ HRP Conjugate, 11 ml	1	pcs	red	until exp. date
5	SUBS TMB Substrate solution, 11 ml	1	pcs	colourless	until exp. date
6	BUF WASH 21X Washing solution concentrate 21x, 22 ml	1	pcs	colourless	Concentrate – until exp. date Diluted washing solution – 1 month at 2...+8 °C or 5 days at RT
7	STOP Stop solution, 11 ml	1	pcs	colourless	until exp. date
8	N003 Plate sealing tape	2	pcs		N/A
9	K206I Instruction Prolactin EIA	1	pcs		N/A
10	K206Q QC data sheet Prolactin EIA	1	pcs		N/A

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 25–250 µl;
- Dry thermostat for +37 °C ±0.1 °C;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

5.3. Storage and stability of the Kit

Store the whole kit at +2...+8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED) to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at +2...+8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

7. TEST PROCEDURE**7.1. Reagent Preparation**

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 21X by 21 dilutions in distilled water.

7.2. Procedural Note:

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1–5 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using (zero calibrator). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
3	Dispense 100 µl of CONJ HRP into the wells.
4	Pipet 50 µl of calibrators CAL 1–6, control samples CONTROL and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
5	Incubate 60 minutes at 37 °C.
6	Prepare washing solution by 21x dilution of washing solution concentrate (BUF WASH 21X) by distilled water. Wash the strips 5 times.
7	Dispense 100 µl of SUBS TMB into the wells
8	Incubate 10–20 minutes at +18...+25 °C
9	Dispense 100 µl of STOP into the wells.
10	Measure OD (optical density) at 450 nm.
11	Set photometer blank on first calibrator
12	Apply point-by-point method for data reduction.

7.5. Handling notes

Calibrators and control sample(s) – only one freezing/thawing cycle is allowed

8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

9. CALCULATION OF RESULTS


9.1. Calculate the mean absorbance values (OD450) for each pair of calibrators and samples.

9.2. Plot a calibration curve on graph paper: OD versus Prolactin concentration.

9.3. Determine the corresponding concentration of Prolactin in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

9.4. Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

Calibrators	Value	Absorbance Units (450 nm)
CAL 1	0 mIU/l	0.06
CAL 2	100 mIU/l	0.33
CAL 3	200 mIU/l	0.61
CAL 4	1000 mIU/l	1.57
CAL 5	2000 mIU/l	2.21

10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for Prolactin. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units, mIU/l	
	Lower limit	Upper limit
Males	60	560
Females		
Pregnancy week:		
1st trimester	-	2000
2nd trimester	-	6000
3rd trimester	-	10000
Menstrual cycle:		
follicular phase	60	600
luteinic phase	120	900
post menopausal	40	550

11. PERFORMANCE CHARACTERISTIC

11.1. Analytical specificity / Cross reactivity

Analyte	Cross-reactivity, % wt/wt
hGH	<0.1
lactogen	<0.1
TSH	<0.1
LH	<0.1
FSH	<0.1

11.2. Analytical sensitivity

Sensitivity of the assay was assessed as being 10 mIU/l.

11.3. Linearity

Linearity was checked by assaying dilution series of 5 samples with different Prolactin concentrations. Linearity percentages obtained ranged within 90 to 110%.

11.4. Recovery

Recovery was estimated by assaying 5 mixed samples with known Prolactin concentrations. The recovery percentages ranged from 90 to 110%.

12. LITERATURE

1. Caufriez A. Menstrual disorders associated with hyperprolactinemia. *Hormone Res.* 1985; 22:209-14.
2. Niall, M.D. et al, «The Chemistry of Growth Hormone and the Lactogenic Hormones»; *Recent Progr. Horm. Res.* 29,471 (1974).
3. Frantz A.G. Physiology in medicine: prolactin. *New Engl. J Med* 1978; 298: 201-7.
4. Kato Y. et al. Regulation of prolactin secretion. In: *The Pituitary Gland*. Ed. H. Imura. New York, Raven Press. 1985; pp.261-278.
5. Thorner, M.O., Edwards, C.R.W., Hanker, J.P., Abraham, G., and Besser, G.M., «The Testes in Normal and Infertile Men», Troen, P. and Nankin, H.R. (eds.), Raven Press, New York, 351-366, (1977).
6. Daughday, W.H., «The Adenohypophysis, Textbook of Endocrinology»; Williams, 6th Ed., Chapter 3, 87-87, (1981).