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1. Oesterling JE Prostate-specific antigen: a critical assessment of the most useful tumor marker fro adenocarcinoma of the prostate. J Urol 1991; 145: 907-23.
2. Oesterling JE et al.. Serum prostate-specific antigen in a community-based population of healthy men: establishment of age-specific reference ranges. JAMA 1993; 270:860.
3. Catalona WJ et al. Measurement of prostate-specific antigen in serum as a screening test for prostate cancer. N Eng J Med 1991, 324: 1156-61.
4. Christensson A et al. Serum prostate-specific antigen complexed to a1-antichymotrypsin as an indicator of prostate cancer. J Urology 1993; 150: 100-5.
5. Milfor Ward A et al. Free/total PSA ratio as an aid to the diagnosis of prostatic carcinoma. Clin Chem 1995; 41:S230.
6. Wu JT. Assay for prostate specific antigen (PSA): problems and possible solutions. J Clin Lab Analysis 1994; 8:51-62.

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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF TOTAL PSA IN HUMAN BLOOD SERUM OR PLASMA

1. INTENDED USE

A solid-phase enzyme immunoassay for the quantitative determination of total PSA in blood serum or plasma.

This kit is designed for measurement of total PSA in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 42 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Prostate specific antigen (PSA) is a serin-like protease with molecular weight ca. 34 kDa and was initially found exclusively in normal prostatic gland as well as in prostatic fluid and seminal plasma. Later it was localized also in breast milk and, according to its enzymological properties, was classified as human prekallikrein 3. In human serum, most of PSA forms complexes with serine protease inhibitor proteins (mostly alpha-1-antichymotrypsin, alpha-2-macroglobulin and antithrypsin). A minor proportion of PSA (free PSA) is circulating outside these complexes.

Elevated serum PSA levels are found in patients with prostatic adenocarcinoma even at early stages of the disease. Values of 1000 ng/ml and even more may be found in patients with profound disease. Clinical value of this parameter is due to possibility of clinical monitoring and prognosis of the disease. Continuous elevation of PSA level is indicative of tumour progression and ineffective therapy. Nevertheless, interpretation of the results obtained should be made in the context of other clinical data. According to data obtained in University of Turku, Finland, the pair of monoclonal antibodies used in present test system (PS2-PS6), recognizes both free and complex-bound forms of PSA with equal affinity (equimolar binding).

Elevations of serum PSA levels are characteristic to prostatic hyperplasia, inflammation and tumours. Serum PSA level can be used for monitoring and treatment control of all diseases involving prostatic tissue, especially prostatic tumours.

Additional information valuable for differential diagnosis between benign and malignant prostate hyperplasia may be obtained by estimation of free PSA/total PSA ratio. In this case, age and case history of patients should be considered: that is, free PSA/total PSA ratio in patients under 60 years is to be estimated if total PSA level is above 4 ng/ml while in males over 60 years when benign prostatic hyperplasia is common this ratio is rational to be estimated when total PSA level is above 7 ng/ml. Besides, it should be kept in mind that significant elevation of total PSA level may be found in patients with prostatitis as well as after massage of prostatic gland and the next day after ejaculation (according to XEMA data, up to 20 ng/ml and 80 ng/ml, respectively).

Please, note, that free PSA/total PSA ratio should be estimated using EIA kits of the same manufacturer. This kit is intended for use with XEMA fPSA EIA, Cat.# K231.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to human total PSA-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Second antibodies – murine monoclonal to human total PSA, labelled with peroxidase enzyme, are then added into the microwells. After washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

4.1. For professional use only.

4.2. This kit is intended for in vitro diagnostic use only.

4.3. INFECTION HAZARD: There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

4.4. Avoid contact with stop solution containing 5,0% H₂SO₄. It may cause skin irritation and burns.

4.5. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

4.6. Do not use the kit beyond the expiration date.

4.7. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

4.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

4.9. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

4.10. Do not mix reagents from different lots.

4.11. Replace caps on reagents immediately. Do not swap caps.

4.12. Do not pipette reagents by mouth.

4.13. Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

4.14. Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

4.15. The Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1 SORB MTP	tpSA EIA strips, 8x12 wells	1	pcs		until exp. date
2 CAL 1-5	polystyrene microwells coated with murine monoclonal to human total PSA human total PSA diluted in tris buffered BSA solution, preservative – 0,01% Bronidox L, 0,01% 2-Methyl-4-isothiazolin-3-one-hydrochloride; also contains red dye	5	pcs	red (C1 – colourless)	2 months
3 CONTROL	dilution of preselected human serum, with high content of total PSA with BSA solution; preservative – 0,01% Bronidox L, 0,01% 2-Methyl-4-isothiazolin-3-one-hydrochloride, colourless	1	pcs	colourless	2 months
4 CONJ HRP	aqueous solution of murine monoclonal to human total PSA coupled with horseradish peroxidase diluted on phosphate buffered solution with casein from bovine milk and detergent (Tween-20), contains 0,1% phenol as preservative and red dye	1	pcs	red	until exp. date
5 SUBS TMB	ready-to-use single-component tetramethylbenzidine (TMB) solution.	1	pcs	colourless	until exp. date
6 BUF WASH 21X	aqueous solution of sodium chloride and detergent (Tween 20), contains proClim300 as a preservative	1	pcs	colourless	Concentrate – until exp. date Diluted washing solution – 1 month at 2...+8°C or 5 days at RT
7 STOP	5,0% vol/vol solution of sulphuric acid	1	pcs	colourless	until exp. date
8 N003	Plate sealing tape	2	pcs		N/A
9 K221I	Instruction tpSA EIA	1	pcs		N/A
10 K221Q	QC data sheet tpSA EIA	1	pcs		N/A

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 25–250 µl;
- Dry thermostat for +37 °C ±0.1 °C;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

5.3. Storage and stability of the Kit

Store the whole kit at +2...+8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED) to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at +2...+8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

7. TEST PROCEDURE**7.1. Reagent Preparation**

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18 to +25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 21X by 21 dilutions in distilled water.

7.2. Procedural Note:

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1-5 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using zero calibrator. Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
3	Dispense 100 µl of CONJ HRP into the wells.
4	Pipet 50 µl of calibrators CAL 1-5, control samples CONTROL and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
5	Incubate 60 minutes at 37 °C.
6	Prepare washing solution by 21x dilution of washing solution concentrate (BUF WASH 21X) by distilled water. Wash the strips 5 times.
7	Dispense 100 µl of SUBS TMB into the wells
8	Incubate 10-20 minutes at +18...+25 °C
9	Dispense 100 µl of STOP into the wells.
10	Measure OD (optical density) at 450 nm.
11	Set photometer blank on first calibrator
12	Apply point-by-point method for data reduction.

7.5. Handling notes

Calibrators and control sample(s) – only one freezing/thawing cycle is allowed

8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

9. CALCULATION OF RESULTS


9.1. Calculate the mean absorbance values (OD450) for each pair of calibrators and samples.

9.2. Plot a calibration curve on graph paper: OD versus total PSA concentration.

9.3. Determine the corresponding concentration of total PSA in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

9.4. Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

Calibrators	Value	Absorbance Units (450 nm)
CAL 1	0 ng/ml	0.04
CAL 2	1.5 ng/ml	0.27
CAL 3	5 ng/ml	0.57
CAL 4	10 ng/ml	0.91
CAL 5	30 ng/ml	2.34

10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for tPSA. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units, ng/ml	
	Lower limit	Upper limit
Males		
<40 yrs	-	4.0
41-60 yrs	-	5.5
>61 yr	-	7.0
Females	-	0.45

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical specificity / Cross reactivity

The pair of monoclonal antibodies used in present test system (PS2-PS6), recognizes both free and complex-bound forms of PSA with equal affinity (equimolar binding)

11.2. Analytical sensitivity. Sensitivity of the assay was assessed as being 0.1 ng/ml.

11.3. Linearity. Linearity was checked by assaying dilution series of 5 samples with different total PSA concentrations. Linearity percentages obtained ranged within 90 to 110%.

11.4. Recovery. Recovery was estimated by assaying 5 mixed samples with known total PSA concentrations. The recovery percentages ranged from 90 to 110%.

12. LITERATURE

- Oesterling JE Prostate-specific antigen: a critical assessment of the most useful tumor marker fro adenocarcinoma of the prostate. J Urol 1991; 145: 907-23.
- Oesterling JE et al.. Serum prostate-specific antigen in a community-based population of healthy men: establishment of age-specific reference ranges. JAMA 1993; 270:860.
- Catalona WJ et al. Measurement of prostate-specific antigen in serum as a screening test for prostate cancer. N Eng J Med 1991, 324: 1156-61.
- Christensson A et al. Serum prostate-specific antigen complexed to a1-antichymotrypsin as an indicator of prostate cancer. J Urology 1993; 150: 100-5.
- Milfor Ward A et al. Free/total PSA ratio as an aid to the diagnosis of prostatic carcinoma. Clin Chem 1995; 41:S230.
- Wu JT. Assay for prostate specific antigen (PSA): problems and possible solutions. J Clin Lab Analysis 1994; 8:51-62.