Instruction for use

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**A SOLID-PHASE ENZYME IMMUNOASSAY
FOR THE QUANTITATIVE DETERMINATION
OF MUC1 (M20) IN HUMAN BLOOD SERUM OR PLASMA**

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CONTENT

| | |
|------------------------------------|----|
| 1. INTENDED USE | 9 |
| 2. SUMMARY AND EXPLANATION | 9 |
| 3. PRINCIPLE OF THE TEST | 10 |
| 4. WARNINGS AND PRECAUTIONS | 10 |
| 5. KIT COMPONENTS | 11 |
| 6. SPECIMEN COLLECTION AND STORAGE | 12 |
| 7. TEST PROCEDURE | 12 |
| 8. QUALITY CONTROL | 14 |
| 9. CALCULATION OF RESULTS | 14 |
| 10. EXPECTED VALUES | 15 |
| 11. PERFORMANCE CHARACTERISTICS | 15 |
| 12. LITERATURE | 15 |

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- McGuire WL, Tandon AK, Allred D, Chamnes GC, Clark GM. How to use prognostic factors in axillary node negative breast cancer patients. J Natl Cancer Inst 1990; 82:1006-7.
- Nicholson S, Richard J, Sainsbury C, et al. Epidermal growth factor receptor (EGFr): results of a 6-year follow up study in operable breast cancer with emphasis on the node-negative subgroup. Br J Cancer 1991; 63:146-50.
- Somerville JE, Clarke LA, Biggart JD. C-erb B-2 overexpression and histological type of in-situ and invasive breast carcinoma. J Clin Pathol 1992;45:16-20.
- Ueronese S, Gambacorta M. Detection of Ki-67 rate in breast cancer. AM J Clin Pathol 1991; 95:30-4.
- Lotnick M, Pavesi F, Scarabelli M. Tumor associated antigens CA15-3 and CA125 in ovarian cancer. Int. J. Biolog Markers 1991; 6:115

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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF MUC1 (M20) IN HUMAN BLOOD SERUM OR PLASMA

1. INTENDED USE

A solid-phase enzyme immunoassay for the quantitative determination of MUC1 (M20) in blood serum or plasma.

This kit is designed for measurement of MUC1 (M20) in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 41 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

MUC1 is a heterogenous glycoprotein with a molecular mass ca. 300–450 kD. Elevation of serum MUC1 is associated with mammary carcinomas. Quantitative determination of MUC1 in serum and plasma is helpful in monitoring of patients with such tumours to estimate the course of the disease, effectiveness of its treatment and to reveal recurrence or metastases. However, MUC1 values obtained should always be interpreted in context of results obtained by other diagnostic and clinical procedures. Besides mammary carcinomas, MUC1 levels in blood may rise in lung tumours, prostate cancer, ovarian carcinomas, gastro-intestinal tumours. Elevation of MUC1 level in blood can be also found in benign tumours of the mammary gland and the ovary, endometriosis, hepatitis, liver cirrhosis and lung fibrosis. Pregnancy and lactation may also cause elevation of MUC1 level in serum.

The M20 kit is based on a principle of competitive immunoassay. The labelled monoclonal antibody X19 recognizing the peptide determinant TRPAPGS binds the natural purified mammary carcinoma-derived antigen on the solid phase. Such test systems (e.g., BR27.29 developed by Biomira) may give a higher sensitivity for tumour specific forms of MUC1 compared to kits based on sandwich principle. However, due to considerable heterogeneity of tumour derived forms of MUC1 antigen, we recommend to use this test in conjunction with two other test systems designed by XEMA for diagnostics and monitoring of mammary carcinoma – M22 (analogous to MCA, Roche) and M12 (analogous to CA15-3 developed by Centocor). All three systems should be used to evaluate MUC1 concentration before and after surgery; the test system showing the most pronounced postsurgery decline should be then used for further monitoring.

3. PRINCIPLE OF THE TEST

This test is based on competition enzyme immunoassay principle. Tested specimen is placed into the microwells coated by purified MUC1 antigen simultaneously with conjugated antibodies. MUC1 antigen from the specimen competes with the immobilized antigen for binding to labelled antibodies. After washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is inversely related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

4.1. For professional use only.

4.2. This kit is intended for in vitro diagnostic use only.

4.3. INFECTION HAZARD: There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

4.4. Avoid contact with stop solution containing 5,0 % H_2SO_4 . It may cause skin irritation and burns.

4.5. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

4.6. Do not use the kit beyond the expiration date.

4.7. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

4.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

4.9. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

4.10. Do not mix reagents from different lots.

4.11. Replace caps on reagents immediately. Do not swap caps.

4.12. Do not pipette reagents by mouth.

4.13. Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

4.14. Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

4.15. The Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5.1. Contents of the Kit

5. KIT COMPONENTS

| Symbol | Description | Qty | Units | Colour code | Stability of opened/diluted components |
|----------------|---|-----|-------|------------------------|---|
| 1 SORB MTP | MUC1 (M20) EIA strips, 8x12 wells | 1 | pcs | | until exp. date |
| 2 CAL 1 – 6 | Calibrator set, 0.8 ml each, zero calibrator C1 – 6 ml. The set contains 6 calibrators; 0; 2,5; 5; 10; 25; 50 Uj/ml | 6 | pcs | blue (C1 – colourless) | 2 months |
| 3 CONTROL | Control serum (0,8 ml) | 1 | pcs | colourless | 2 months |
| 4 CONJ HRP | Conjugate, 6 ml | 1 | pcs | red | until exp. date |
| 5 SUBS TMB | Substrate solution, 11 ml | 1 | pcs | colourless | until exp. date |
| 6 BUF WASH 21X | Washing solution concentrate 21x, 22 ml | 1 | pcs | colourless | Concentrate – until exp. date Diluted washing solution – 1 month at 2-8 °C or 5 days at RT |
| 7 STOP | Stop solution, 11 ml | 1 | pcs | colourless | until exp. date |
| 8 N003 | Plate sealing tape | 2 | pcs | | N/A |
| 9 K228I | Instruction MUC1 (M20)EIA | 1 | pcs | | N/A |
| 10 K228Q | QC data sheet MUC1 (M20)EIA | 1 | pcs | | N/A |

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 25–250 µl;
- Dry thermostat for 37 °C ±0.1 °C;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

5.3. Storage and stability of the Kit

Store the whole kit at 2 to 8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells **TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED)** to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at +2...+8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

7. TEST PROCEDURE**7.1.** Reagent Preparation

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 21X by 21 dilutions in distilled water.

7.2. Procedural Note:

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Assay procedure

| | |
|----|--|
| 1 | Put the desired number of microstrips into the frame; allocate 14 wells for the calibrators CAL 1–6 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS. |
| 2 | If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using zero calibrator. Use of other buffers or reagents for sample dilution may lead to incorrect measurement. |
| 3 | Pipet 50 µl of calibrators CAL 1–6, control samples CONTROL and unknown samples into the wells. |
| 4 | Dispense 50 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape (included into the kit). |
| 5 | Incubate 60 minutes at 37 °C. |
| 6 | Prepare washing solution by 21x dilution of washing solution concentrate (BUF WASH 21X) by distilled water. Wash the strips 5 times. |
| 7 | Dispense 100 µl of SUBS TMB into the wells |
| 8 | Incubate 10-20 minutes at 18-25 °C |
| 9 | Dispense 100 µl of STOP into the wells. |
| 10 | Measure OD (optical density) at 450 nm. |
| 11 | Set photometer blank on air |
| 12 | Apply lin-log method for data reduction. |

7.5. Handling notes

Calibrators and control sample(s) – only one freezing/thawing cycle is allowed

8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer’s instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

9. CALCULATION OF RESULTS


9.1. Calculate the mean absorbance values (OD450) for each pair of calibrators and samples.

9.2. Plot a calibration curve on graph paper: OD versus MUC1 (M20) concentration.

9.3. Determine the corresponding concentration of MUC1 (M20) in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

9.4. Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

| Calibrators | Value | Absorbance Units (450 nm) |
|-------------|----------|---------------------------|
| CAL 1 | 0 U/ml | 2.28 |
| CAL 2 | 2.5 U/ml | 1.40 |
| CAL 3 | 5 U/ml | 1.12 |
| CAL 4 | 10 U/ml | 0.84 |
| CAL 5 | 25 U/ml | 0.55 |
| CAL 6 | 50 U/ml | 0.36 |

10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for MUC1 (M20). Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

| Sex, age | Units, U/ml | |
|-----------------|-------------|-------------|
| | Lower limit | Upper limit |
| Males | - | 10.0 |
| Females | - | 10.0 |
| Pregnancy week: | | |
| 1st trimester | - | 12 |
| 2nd trimester | 2.0 | 14 |
| 3rd trimester | 2.0 | 18 |
| Lactation | - | 16 |

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical specificity / Cross reactivity.

Sensitivity of the assay was assessed as being 0.5 U/ml.

11.2. Linearity.

Linearity was checked by assaying dilution series of 5 samples with different MUC1 (M20) concentrations. Linearity percentages obtained ranged within 90 to 110 %.

11.3. Recovery.

Recovery was estimated by assaying 5 mixed samples with known MUC1 (M20) concentrations. The recovery percentages ranged from 90 to 110 %.

12. LITERATURE

- McGuire WL, Tandon AK, Allred D, Chamnes GC, Clark GM. How to use prognostic factors in axillary node negative breast cancer patients. J Natl Cancer Inst 1990; 82:1006-7.
- Nicholson S, Richard J, Sainsbury C, et al. Epidermal growth factor receptor (EGFr): results of a 6-year follow up study in operable breast cancer with emphasis on the node-negative subgroup. Br J Cancer 1991; 63:146-50.
- Somerville JE, Clarke LA, Biggart JD. C-erb B-2 overexpression and histological type of in-situ and invasive breast carcinoma. J Clin Pathol 1992;45:16-20.
- Ueronese S, Gambacorta M. Detection of Ki-67 rate in breast cancer. AM J Clin Pathol 1991; 95:30-4.
- Lotnicker M, Pavesi F, Scarabelli M. Tumor associated antigens CA15-3 and CA125 in ovarian cancer. Int. J. Biolog Markers 1991; 6:115

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