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1. Steven Black, Irving Kushner, and David Samols – C-reactive Protein. J. Biol. Chem., Nov 2004; 279: 48487 – 48490.
2. V B Patel, M A Robbins, and E J Topol – C-reactive protein: a 'golden marker' for inflammation and coronary artery disease. Cleveland Clinic Journal of Medicine, Jun 2001; 68: 521 – 524.
3. Nader Rifai and Paul M. Ridker – High-Sensitivity C-Reactive Protein: A Novel and Promising Marker of Coronary Heart Disease. Clin. Chem., Mar 2001; 47: 403 – 411.

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Instruction for use

**A SOLID-PHASE ENZYME IMMUNOASSAY
FOR THE QUANTITATIVE DETERMINATION
OF C-REACTIVE PROTEIN
IN HUMAN BLOOD SERUM OR PLASMA**

1. INTENDED USE

A solid-phase enzyme immunoassay for the quantitative determination of C-reactive protein in blood serum or plasma.

This kit is designed for measurement of C-reactive protein in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 41 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

C-reactive protein (CRP) is a member of pentraxin family with the MW ca. 100 kDa. In evolution of vertebrates, CRP was one of the first developed factors of inborn humoral immunity. CRP efficiently binds both endogenous (phosphoethanolamine, phosphorylcholine, histones, fibronectin, laminin) and exogenous (cell walls of Gram positive bacteria) ligands. The binding of these ligands as well as the debris of apoptotic and necrotic cells by CRP facilitates the phagocytosis. CRP is capable to activate classic pathway of complement activation and phagocytic activity of macrophages. Dramatic (up to 1000 fold) elevation of serum CRP level is highly sensitive but nonspecific marker of inflammation induced by interleukin 6.

Recently, there was observed a correlation between moderately elevated serum CRP (3-10 mg/l) and risk of ischemic heart disease. Elevated serum CRP is also observed in smokers and post-menopausal women receiving substitution hormonal therapy.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to human CRP-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Unbound material is removed by washing procedure. Second antibodies – murine monoclonal to human CRP, labelled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

4.1. For professional use only.

4.2. This kit is intended for in vitro diagnostic use only.

4.3. INFECTION HAZARD: There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

4.4. Avoid contact with stop solution containing 5,0% H_2SO_4 . It may cause skin irritation and burns.

4.5. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

4.6. Do not use the kit beyond the expiration date.

4.7. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

4.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

4.9. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

4.10. Do not mix reagents from different lots.

4.11. Replace caps on reagents immediately. Do not swap caps.

4.12. Do not pipette reagents by mouth.

4.13. Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

4.14. Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

4.15. The Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	CRP EIA strips, 8x12 wells	1	pcs		until exp.date
2	polystyrene microwells coated with murine monoclonal to human CRP	6	pcs	bright red (C1 - colourless)	2 months
3	human C-reactive protein diluted in tris buffered BSA solution, preservative - 0,01% Bronidox L, 0,01% 2-Methyl-4-isothiazolin-3-one-hydrochloride; also contains bright red dye	1	pcs	colourless	2 months
4	dilution of preselcted human serum, with high content of C-reactive protein with BSA solution; preservative - 0,01% Bronidox L, 0,01% 2-Methyl-4-isothiazolin-3-one-hydrochloride, colourless	1	pcs	bright red	until exp.date
5	aqueous solution of murine monoclonal to human CRP coupled with horseradish peroxidase diluted on phosphate buffered solution with casein from bovine milk and detergent (Tween-20), contains 0,1% phenol as preservative and bright red dye	1	pcs	red	until exp.date
6	phosphate buffered saline with casein from bovine milk and detergent (Tween-20), contains 0,1% phenol as preservative; contains red dye	1	pcs	blue	until exp.date
7	phosphate buffered saline with casein from bovine milk and detergent (Tween-20), contains 0,1% phenol as preservative and blue dye	1	pcs	colourless	until exp.date
8	ready-to-use single-component tetramethylbenzidine (TMB) solution.	1	pcs	colourless	Concentrate - until exp.date Diluted washing solution - 1 month at 2...+8 °C or 5 days at RT
9	aqueous solution of sodium chloride and detergent (Tween 20), contains proClin300 as a preservative	1	pcs	colourless	until exp.date
10	5,0% vol/vol solution of sulphuric acid	2	pcs		N/A
11	Stop solution, 11 ml	1	pcs		N/A
12	Plate sealing tape	1	pcs		N/A
	Instruction CRP EIA	1	pcs		N/A
	QC data sheet CRP EIA	1	pcs		N/A

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 μl , is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 5–250 μl ;
- Dry thermostat for 37 $^{\circ}\text{C} \pm 0.1$ $^{\circ}\text{C}$
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

5.3. Storage and stability of the Kit

Store the whole kit at +2...+8 $^{\circ}\text{C}$ upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED) to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at +2...+8 $^{\circ}\text{C}$ before testing. For a longer storage, the specimens should be frozen at -20 $^{\circ}\text{C}$ or lower. Repeated freezing/thawing should be avoided.

7. TEST PROCEDURE**7.1.** Reagent Preparation

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25 $^{\circ}\text{C}$) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 21X by 21 dilutions in distilled water.

7.2. Procedural Note:

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 14 wells for the calibrators CAL 1–6 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Dilute all samples using buffer DIL SPE (EIA buffer) 101 fold. Do not dilute control sample and calibrators.
3	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using DIL SPE (EIA buffer). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
4	Pipet 100 µl of red EIA buffer into each well.
5	Pipet 25 µl of calibrators CAL 1–6, control samples CONTROL and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
6	Incubate 30 minutes at 37 °C .
7	Prepare washing solution by 21x dilution of washing solution concentrate BUF WASH 21X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times
8	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
9	Incubate 30 minutes at 37 °C .
10	Wash the strips 5 times.
11	Dispense 100 µl of SUBS TMB into the wells
12	Incubate 10–20 minutes at +18...+25 °C
13	Dispense 100 µl of STOP into the wells.
14	Measure OD (optical density) at 450 nm.
15	Set photometer blank on first calibrator
16	Apply point-by-point method for data reduction.

7.5. Sample processing

Material type	Notes on material collection, storage and handling	Sample dilution example	red EIA buffer into the well, µl	Sample into the well, µl	Calculation factor
blood serum or plasma	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	5 µl of sample + 500 µl of diluent	100	25	1

8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

9. CALCULATION OF RESULTS


9.1. Calculate the mean absorbance values (OD450) for each pair of calibrators and samples.

9.2. Plot a calibration curve on graph paper: OD versus C-reactive protein concentration.

9.3. Determine the corresponding concentration of C-reactive protein in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

9.4. Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

Calibrators	Value	Absorbance Units (450 nm)
CAL 1	0 mg/l	0.05
CAL 2	0.25 mg/l	0.11
CAL 3	1 mg/l	0.32
CAL 4	2.5 mg/l	0.74
CAL 5	10 mg/l	1.97
CAL 6	25 mg/l	2.74

10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for CRP. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units, mg/l	
	Lower limit	Upper limit
Healthy donors	-	5.0

11. PERFORMANCE CHARACTERISTICS**11.1. Analytical sensitivity**

Sensitivity of the assay was assessed as being 0.05 mg/l.

11.2. Linearity

Linearity was checked by assaying dilution series of 5 samples with different C-reactive protein concentrations. Linearity percentages obtained ranged within 90 to 110%.

11.3. Recovery

Recovery was estimated by assaying 5 mixed samples with known C-reactive protein concentrations. The recovery percentages ranged from 90 to 110%.

12. LITERATURE

1. Steven Black, Irving Kushner, and David Samols – C-reactive Protein. *J. Biol. Chem.*, Nov 2004; 279: 48487 – 48490.
2. V B Patel, M A Robbins, and E J Topol – C-reactive protein: a 'golden marker' for inflammation and coronary artery disease. *Cleveland Clinic Journal of Medicine*, Jun 2001; 68: 521 – 524.
3. Nader Rifai and Paul M. Ridker – High-Sensitivity C-Reactive Protein: A Novel and Promising Marker of Coronary Heart Disease. *Clin. Chem.*, Mar 2001; 47: 403–411.