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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF TOTAL IgG IN HUMAN BIOLOGICAL FLUIDS

1. INTENDED USE

A solid-phase enzyme immunoassay for the quantitative determination of total IgG in biological fluids.

This kit is designed for measurement of total IgG in biological fluids. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 42 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Immunoglobulin G (IgG) is the main part of serum γ - globulin fraction. IgG is secreted during secondary immune response and plays a key role in humoral immunity. Decrease of serum IgG concentration below 5 g/l is a marker of severe life-threatening immunodeficiency. Determination of serum IgG concentration and IgG/IgA/IgM ratios can be used for monitoring of humoral immune status. Marked elevation of serum IgG may be observed in chronic inflammation, autoimmune diseases and myeloma.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to human total IgG-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Unbound material is removed by washing procedure. Second antibodies – murine monoclonal to human total IgG, labelled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

4.1. For professional use only.

4.2. This kit is intended for in vitro diagnostic use only.

4.3. INFECTION HAZARD: There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

4.4. Avoid contact with stop solution containing 5,0% H_2SO_4 . It may cause skin irritation and burns.

4.5. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

4.6. Do not use the kit beyond the expiration date.

4.7. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

4.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

4.9. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

4.10. Do not mix reagents from different lots.

4.11. Replace caps on reagents immediately. Do not swap caps.

4.12. Do not pipette reagents by mouth.

4.13. Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

4.14. Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

4.15. The Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	SORB MTP total IgG EIA strips, 8x12 wells	1	pcs		until exp.date
2	CAL 1-5 Calibrator set, 1 ml each. The set contains 5 calibrators: 0; 1; 5; 10; 25 g/l	5	pcs	blue (Cl - colourless)	2 months
3	CONTROL Control serum (1 ml)	1	pcs	colourless	2 months
4	CONJ HRP Conjugate, 11 ml	1	pcs	red	until exp.date
5	DIL SPE EIA sample buffer 100 ml	1	pcs	red	until exp.date
6	SUBS TMB Substrate solution, 11 ml	1	pcs	red	until exp.date
7	BUF WASH 21X Washing solution concentrate 21x, 22 ml	1	pcs	colourless	Concentrate - until exp.date Diluted washing solution - 1 month at 2...+8 °C or 5 days at RT
8	STOP Stop solution, 11 ml	1	pcs	colourless	until exp.date
9	N003 Plate sealing tape	2	pcs		N/A
10	K2711 Instruction total IgG EIA	1	pcs		N/A
11	K271Q QC data sheet total IgG EIA	1	pcs		N/A

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 μl , is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 10–250 μl ;
- Dry thermostat for 37 $^{\circ}\text{C} \pm 0.1$ $^{\circ}\text{C}$
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

5.3. Storage and stability of the Kit

Store the whole kit at +2...+8 $^{\circ}\text{C}$ upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED) to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at +2...+8 $^{\circ}\text{C}$ before testing.

7. TEST PROCEDURE**7.1. Reagent Preparation**

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25 $^{\circ}\text{C}$) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 21X by 21 dilutions in distilled water.

7.2. Procedural Note:

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1-5 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Dilute samples using buffer DIL SPE (EIA sample buffer) 5000 fold. See table M for dilution modes and factors for different types of analyzed material. Do not dilute control sample and calibrators.
3	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using DIL SPE (EIA sample buffer). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
4	Pipet 100 µl of calibrators CAL 1-5 and control samples CONTROL into allocated wells. For testing of blood serum or plasma pipet 100 µl of the unknown diluted sample (DILUTION 2) into the allocated wells. See table M for the volumes of other materials. Pipetting should be made within 3 minutes, to ensure an uniform incubation time for all samples. Carefully mix the contents of the wells by short horizontal rotating of the plate for 5-7 seconds and cover the wells by plate adhesive tape (included into the kit).
5	Incubate 30 minutes at +37 °C.
6	Prepare washing solution by 21x dilution of washing solution concentrate BUF WASH 21X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
7	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
8	Incubate 30 minutes at +37 °C.
9	Wash the strips 5 times.
10	Dispense 100 µl of SUBS TMB into the wells
11	Incubate 10-20 minutes at +18...+25 °C
12	Dispense 100 µl of STOP into the wells.
13	Measure OD (optical density) at 450 nm.
14	Set photometer blank on first calibrator
15	Apply point-by-point method for data reduction. Use Calculation factor listed in table M to calculate analyte concentration in different material types.

7.5. Sample processing

Material type	Notes on material collection, storage and handling	Sample dilution example	EIA sample buffer into the well, μl	Sample into the well, μl	Calculation factor
blood serum or plasma	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	10 μl of sample + 990 μl of diluent = DILUTION 1. 10 μl of DILUTION 1 + 490 μl of diluent = DILUTION 2	0	100	1
saliva	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		90	10	0.002
urine	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		50	50	0.0004
cerebrospinal fluid	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	10 μl of sample + 500 μl of diluent	0	100	0.01

8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

9. CALCULATION OF RESULTS

9.1.1. Calculate the mean absorbance values (OD450) for each pair of calibrators and samples.


9.2. Plot a calibration curve on graph paper: OD versus total IgG concentration.

9.3. Determine the corresponding concentration of total IgG in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is

recommended due to non-linear shape of curve.

9.4. Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

Calibrators	Value	Absorbance Units (450 nm)
CAL 1	0 g/l	0.08
CAL 2	1 g/l	0.38
CAL 3	5 g/l	1.15
CAL 4	10 g/l	1.78
CAL 5	25 g/l	2.30

10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for total IgG. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units, g/l	
	Lower limit	Upper limit
newborn	7.0	15
1-3 month	2.7	8.0
4-6 month	1.8	8.5
7-12 month	3.5	12
1-6 yrs	6.5	18
7-11 yrs	8.5	15
> 11 yrs	9.0	20

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical specificity / Cross reactivity

Analyte	Cross-reactivity, % wt/wt
IgA	<0.1
IgM	<0.1
IgE	<0.1

11.2. Analytical sensitivity. Sensitivity of the assay was assessed as being 0.12 g/l.

11.3. Linearity. Linearity was checked by assaying dilution series of 5 samples with different total IgG concentrations. Linearity percentages obtained ranged within 90 to 110%.

11.4. Recovery. Recovery was estimated by assaying 5 mixed samples with known total IgG concentrations. The recovery percentages ranged from 90 to 110%.

12. LITERATURE

- RG Hamilton – Human IgG subclass measurements in the clinical laboratory. Clin. Chem., Oct 1987; 33: 1707 – 1725.
- V. A. Semenova, E. Steward-Clark, K. L. Stamey, T. H. Taylor, Jr., D. S. Schmidt, S. K. Martin, N. Marano, and C. P. Quinn – Mass Value Assignment of Total and Subclass Immunoglobulin G in a Human Standard Anthrax Reference Serum. Clin. Diagn. Lab. Immunol., Sep 2004; 11: 919 – 923.