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1. Wisdom GB, Enzyme-Immunoassay, Clin. Chem.22: 1243-1255,1976.
2. Wolters G, Kuipers L, Kacaki J, and Schuurs A, Solid-phase enzyme-immunoassay for detection of hepatitis B surface antigen, J. Clin. Pathol. 29:873-879
3. Wei R, Lmogjt GJ, Zimmerman DH, and Bond HE, Solid-Phase Enzyme immunoassay for Hepatitis B Surface Antigen, Clin. Chem., 23:813-815, 1977.
4. David GS, Present W, Martinis J, Wang R, Bartholomew R, Desmond W, and Sevier ED, Monoclonal antibodies in the detection of hepatitis infection, Med. Lab. Sci. 38:341-348, 1981.
5. Goodall AH, Miescher G, Meek FM, Janossy G, Thomas HC, Monoclonal antibodies in a solid-phase radiometric assay for HBsAg Med. Lab. Sci. 38:349-354, 1981
6. Kennedy RC, Ionscu-Matiu I, Alder-Storthz K, Henkel RD, Sanchez Y, Dreesman GR, Characterization of Anti-Hepatitis B Surface Antigen Monoclonal Antibodies, Intervirology. 19:176-180, 1983.
7. Shih JW-K, Cote PJ, Dapolito GM, and Gerin JL, Production of monoclonal antibody against Hepatitis B surface antigen(HBsAg) by somatic cell hybrids, J Virol. Meth. 1:257-273,1980.
8. Wands JR, Zurawski VR, High Affinity Monoclonal Antibodies to Hepatitis B Surface Antigen(HBsAg) Produced by Somatic Cell Hybrids, Gastroenterology 80:225-232,1981
9. S. Department of Health and Human Services. Biosafety in Microbiological and biomedical laboratories. HHS Publication(NIH) 88-8395. Washington: U.S. Government Printing Office, May 1988.
10. World Health Organization. Laboratory Biosafety manual. Geneva. World Health Organization, 1983.
11. National Committee for Clinical Laboratory Standards. Protection of laboratory workers from infectious disease transmitted by blood, body fluids, and tissue: Tentative guideline. NCCLS Document M29-T. Villanova, PA.: NCCLS,1989.
12. Centers for Disease Control. Recommendation for prevention of HIV transmission in health care setting. MMWR 36, Supplement No. 2S, 1987.
13. Schulster, L. M., Hollinger, F. B., Dreesman. G. R., and Melnick, J. L., Immunological and biophysical alteration of Hepatitis B virus antigens by sodium hypochlorite disaffection. Appl. And Envir. Microbiol., 42:762-767, 1981.

**InTec PRODUCTS, INC.**

Вироблено для New Vision Solutions Limited