PRESTIGE 24i LQ LIPASE

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1. Tietz NW et al. Lipase in serum-the elusive enzyme: An overview. Clin Chem 1993;39:746-756.
2. Steinberg WM, Goldstein SS, Davies ND et al. Diagnostic assays in acute pancreatitis. (Review). Ann Intern Med 1985; 102:576-580.
3. Leybold A, Junge W. Importance of colipase for the measurement of serum lipase activity. Adv clin Enzymol 1986;4:60-67.
4. Alan H. B. Wu, Tietz Clinical Guide to Laboratory Tests, W.B. Saunders Company, 4th edition, 676 (2006).

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Item name	26	LIPA				
Auto Rerun SW						
ON						
Auto Rerun Range (Result)						
	ON	ON				
	Lower	Higher				
Serum	7	600				
Urine						
Plasma						
CSF						
Dialysis						
Other						
Auto Rerun Condition (Absorbance)						
Absorbance Range						
	Lower	OFF				
	Higher	OFF				
Prozone Range						
OFF						

АДАПТАЦИЯ для Prestige 24i Premium, Biolis 24i Premium и Sapphire 400 Premium

Item No.	26	Item Name	LIPA				
Data information							
Units	U/l						
Decimals	1						
Calibration							
Type	Linear2						
Std sample conc.							
Blank	0	#1	*	#2	*		
#3		#4		#5			
#6							
Analysis							
Type	RATE method						
Main Wave Length	570 nm						
Sub Wave Length	700 nm						
Method	Colorimetric						
Correlation							
	Slope	X+		Intercept			
Y=	1			0			

Item No.	26	Item Name	LIPA				
Aspiration							
Kind	Double						
Data Process							
Read	Start	End					
Main	35	42					
Sub							
Abs.Limit							
	Low	High					
	-3	~		3			
Blank value							
Water Blank							
Reaction Monitor							
0 Level Point	1						
Span	3						
Third mixing							
OFF							
Correction value							
Blank correction							
End Point Limit	2						
Linear Check (%)	80						
Prozone Check							
	Start	End	Limit (%)				
First							
Second			Low				

Item No.	26	Item Name	LIPA				
Normal Range							
	Male			Female			
	Low	High	Low	High	Low	High	
Serum	13	60	13	60			
Urine							
Plasma	13	60	13	60			
CSF							
Dialysis							
Other							
Panic Range							
	Male			Female			
	Low	High	Low	High	Low	High	
Serum							
Urine							
Plasma							
CSF							
Dialysis							
Other							

Item No.	26	Item Name	LIPA				
Auto Rerun SW							
ON							
Auto Rerun Range (Conc.)							
	First Dil	Low			High		
	Re	Value	Dil	Re	Value	Dil	
Serum		6			1000		
Urine							
Plasma							
CSF							
Dialysis							
Other							
Auto Rerun Condition (Absorbance)							
Lower	OFF						
Higher	OFF						
Auto Rerun Condition (Prozone)							
OFF							
Dilution							
100:Dil2							

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PZ CORMAY S.A.