PRESTIGE 24i HDL DIRECT

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1. National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
2. Shih WJ, Bachorik PS, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3:351 – 364.
3. Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III), NIH Publication No. 02-5215, September 2002.
4. Badimon J.J., Badimon L., Fuester V.: Journal of Clinical Investigation, 85:1234-41 (1990).
5. Castelli W.P. et al.: Circulation, 55; 767 (1977).
6. Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 564 (2006).

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Item name	6	HDL-D
Auto Rerun SW		
ON		
Auto Rerun Range (Result)		
	ON	ON
	Lower	Higher
Serum	1.6	465
Urine		
Plasma		
CSF		
Dialysis		
Other		
Auto Rerun Condition (Absorbance)		
Absorbance Range		
Lower	OFF	
Higher	OFF	
Prozone Range		
	OFF	

АДАПТАЦИЯ для Prestige 24i Premium, Biolis 24i Premium и Sapphire 400 Premium

Item No.	6	Item Name	HDL-D	Optical
Data information				
Units	mg/dl			
Decimals	1			
Analysis				
Type	END method			
Main Wave Length	600 nm			
Sub Wave Length	700 nm			
Method	Direct			
Correlation				
	Slope	Intercept		
Y=	1	X+ 0		
Calibration				
Type	Linear1			
Std sample conc.				
Blank	0	#1	*	#2
#3		#4		#5
#6				

Item No.	6	Item Name	HDL-D	Optical
Aspiration				
Kind	Double			
Vol.				
Kind	Vol.	Add	Units	
Sample	3	5	µl	
Reagent 1	270	10	µl	
Reagent 2	90	10	µl	
Data Process				
Read	Start	End		
Main	51	52		
Sub	29	31		
Abs.Limit				
Low	High			
-0.2	~			2
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
OFF				
Correction value				
Blank correction				
End Point Limit				
2				
Linear Check (%)				
Prozone Check				
	Start	End	Limit (%)	
First				
Second			Low	

Item No.	6	Item Name	HDL-D	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	40	60	40	60
Urine				
Plasma	40	60	40	60
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	6	Item Name	HDL-D	Optical
Auto Rerun SW				
ON				
Auto Rerun Range (Conc.)				
	First	Low		High
	Dil	Re	Dil	Re
Serum		Value	Dil	Value
		1.3		280
Urine				
Plasma				
CSF				
Dialysis				
Other				
Auto Rerun Condition (Absorbance)				
Lower	OFF			
Higher	OFF			
Auto Rerun Condition (Prozone)				
	OFF			
Dilution				
	100:Di12			

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PZ CORMAY S.A.