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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF TOTAL IgE IN HUMAN BLOOD SERUM OR PLASMA

1. INTENDED USE

A solid-phase enzyme immunoassay for the quantitative determination of total IgE in blood serum or plasma.

This kit is designed for measurement of total IgE in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 42 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Total immunoglobulin E (IgE) serum level is widely reported as the laboratory marker of atopic diseases such as atopic asthma, atopic dermatitis, and pollenosis. Separately, high levels of total serum IgE are characteristic for parasitic infestations and some other clinical disorders including superficial and systemic mycosis. Decreased levels of IgE are found in cases of hypogammaglobulinemia, autoimmune diseases, ulcerative colitis, and primary biliary cirrhosis.

In allergic patients, serum total IgE level in general corresponds to the severity of the allergic disease and may be therefore used for monitoring of all kinds of anti-allergic therapy or allergen elimination.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific mAb to IgE epsilon. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Unbound material is removed by washing procedure. Second antibodies – rabbit polyclonal antibodies to IgE, labelled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

4.1. For professional use only.

4.2. This kit is intended for in vitro diagnostic use only.

4.3. INFECTION HAZARD: There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

4.4. Avoid contact with stop solution containing 5,0% H_2SO_4 . It may cause skin irritation and burns.

4.5. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

4.6. Do not use the kit beyond the expiration date.

4.7. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

4.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

4.9. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

4.10. Do not mix reagents from different lots.

4.11. Replace caps on reagents immediately. Do not swap caps.

4.12. Do not pipette reagents by mouth.

4.13. Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

4.14. Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

4.15. The Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1 SORB MTP	Total IgE EIA strips, 8x12 wells	1	pcs		until exp.date
2 CAL 1-5	Calibrator set, 0.8 ml each. The set contains 5 calibrators: 0; 50; 200; 500; 1000 IU/ml	5	pcs	red (C1 – colourless)	2 months
3 CONTROL	Control serum (0.8 ml)	1	pcs	colourless	2 months
4 CONJ HRP	Conjugate, 11 ml	1	pcs	red	until exp.date
5 DIL	EIA buffer, 11 ml	1	pcs	blue	until exp.date
6 SUBS TMB	Substrate solution, 11 ml	1	pcs	colourless	until exp.date
7 BUF WASH 21X	Washing solution concentrate 21x, 22 ml	1	pcs	colourless	Concentrate – until exp.date Diluted washing solution – 1 month at 2...+8°C or 5 days at RT
8 STOP	Stop solution, 11 ml	1	pcs	colourless	until exp.date
9 N003	Plate sealing tape	2	pcs		N/A
10 K200I	Instruction Total IgE EIA	1	pcs		N/A
11 K200Q	QC data sheet Total IgE EIA	1	pcs		N/A

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 25–250 µl;
- Dry thermostat for 37 °C ±0,1 °C;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0.

5.3. Storage and stability of the Kit

Store the whole kit at 2 to 8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells **TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED)** to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at 2...+8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

7. TEST PROCEDURE**7.1. Reagent Preparation**

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 21X by 21 dilutions in distilled water.

7.2. Procedural Note:

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Alternative units:

1 IU/ml = 2.15 pg/ml

7.5. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1-5 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS. NOTE: the calibrator/control and unknown sample wells are filled differently.
2	Pipet 50 µl of EIA buffer into the wells allocated for calibrator and controls. Pipet 50 µl of EIA buffer into the wells allocated for samples.
3	Pipet 50 µl of calibrators CAL 1-5, control sample CONTROL and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
4	Incubate 30 minutes at 37 °C.
5	Prepare washing solution by 21x dilution of washing solution concentrate BUF WASH 21X by distilled water; Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
6	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape
7	Incubate 30 minutes at 37 °C.
8	Wash the strips 5 times.
9	Dispense 100 µl of SUBS TMB into the wells
10	Incubate 10-20 minutes at +18...+25 °C.
11	Dispense 100 µl of STOP into the wells.
12	Measure OD (optical density) at 450 nm.
13	Set photometer blank on first calibrator
14	Apply point-by-point method for data reduction.

7.6. Handling notes

Incubation at 37 °C may be replaced by incubation at +20...+25 °C with continuous shaking

8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.


The test must be performed exactly as per the manufacturer’s instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

9. CALCULATION OF RESULTS

1. Calculate the mean absorbance values (OD450) for each pair of calibrators and samples.
2. Plot a calibration curve on graph paper: OD versus total IgE concentration.
3. Determine the corresponding concentration of total IgE in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.
4. Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

Calibrators	Value	Absorbance Units (450 nm)
CAL 1	0 IU/ml	0.08
CAL 2	50 IU/ml	0.36
CAL 3	200 IU/ml	0.75
CAL 4	500 IU/ml	1.32
CAL 5	1000 IU/ml	2.08

10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for Total IgE. Based on data obtained by XEMA, the following normal range is recommended (see below).

Sex, age	Units, IU/ml		Units alternative, pg/ml	
	Lower limit	Upper limit	Lower limit	Upper limit
<6 months	-	12	-	25.8
6-12 months	-	30	-	64.5
1-3 years	-	45	-	96.8
4-6 years	-	70	-	150.5
7-9 years	-	90	-	193.5
10-15 years	-	120	-	258
>15 years	-	130	-	279.5

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical specificity / Cross reactivity

Analyte	Cross-reactivity, % wt/wt
IgG	<0.1
IgM	<0.1
IgA	<0.1

11.2. Analytical sensitivity.

Sensitivity of the assay was assessed as being 6 IU/ml.

11.3. Linearity.

Linearity was checked by assaying dilution series of 5 samples with different total IgE concentrations. Linearity percentages obtained ranged within 90 to 110%.

11.4. Recovery.

Recovery was estimated by assaying 5 mixed samples with known total IgE concentrations. The recovery percentages ranged from 90 to 110%.

12. LITERATURE

1. Zetterstrom and Hohansson S.G.O. Allergy 1981; 36:537.
2. Buckley R. H. Immunopharmacology of Allergic Disease 1979; 117.
3. Michel f. B., Bousquet J. and Greillier P. J. Allergy Clin. Immunol. 1980; 64:422.
4. Ishizaka T. Ann Allergy 1982; 48: 313.
5. Kulczynski A. Jr. J. Allergy Clin. Immunol. 1981; 68:5.