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*Instruction for use*

## **A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF CORTISOL IN HUMAN BLOOD SERUM OR PLASMA**

### **1. INTENDED USE**

A solid-phase enzyme immunoassay for the quantitative determination of cortisol in blood serum or plasma.

This kit is designed for measurement of cortisol in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 41 unknown samples in duplicates.

### **2. SUMMARY AND EXPLANATION**

Cortisol is a glucocorticoid with a MW of 362.5 Dalton. Cortisol is the major hormone secreted by adrenals. In blood, cortisol is found mostly in a bound form, transcortin being the carrier. Cortisol secretion undergoes circadian rhythms with maximal (up to 700 nmol/l) concentrations seen in the morning (6–9 am) and minimal (up to 55 nmol/l) – in the midnight.

During pregnancy, Cortisol blood level is continuously increasing by up to 5-fold of initial concentration before delivery, its circadian rhythm being altered. Cortisol plays an important role in development of alveolar epithelium and surfactant secretion, this being of major importance for the first inhale of a newborn.

Elevated Cortisol concentrations in blood are found in secreting tumours of adrenals, in virilizing hyperplasia of adrenals, in Cushing syndrome, in ACTH-producing tumours, during surgical stress, in cardiac insufficiency, diabetes, burns, pains, during pregnancy, during estrogen therapy, etc. Cortisol blood level may be increased by intake of ACTH, Cortisol, alcohol, nicotine, oral contraceptives.

Decreased Cortisol levels are found in Addison syndrome, adrenogenital syndrome, hypopituitarism. Some drugs may decrease Cortisol level in blood, such as: L-DOPA, dexamethasone, etc. Decreased Cortisol level during pregnancy may indicate anencephaly of the fetus.

### 3. PRINCIPLE OF THE TEST

This test is based on competition enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to cortisol-antibodies simultaneously with conjugated Cortisol-peroxidase. Cortisol from the specimen competes with the conjugated Cortisol for coating antibodies. After washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is inversely related to the quantity of the measured analyte in the specimen.

### 4. WARNINGS AND PRECAUTIONS

**4.1.** For professional use only.

**4.2.** This kit is intended for in vitro diagnostic use only.

**4.3. INFECTION HAZARD:** There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

**4.4.** Avoid contact with stop solution containing 5,0% H<sub>2</sub>SO<sub>4</sub>. It may cause skin irritation and burns.

**4.5.** Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

**4.6.** Do not use the kit beyond the expiration date.

**4.7.** All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

**4.8.** Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

**4.9.** Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

**4.10.** Do not mix reagents from different lots.

**4.11.** Replace caps on reagents immediately. Do not swap caps.

**4.12.** Do not pipette reagents by mouth.

**4.13.** Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

**4.14.** Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

**4.15.** The Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

## 5.1. Contents of the Kit

## 5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	SORB MTP Cortisol EIA strips, 8x12 wells	1	pcs		until exp. date
2	CAL 1-6 Calibrator set, 0.8 ml each. The set contains 6 calibrators: 0; 40; 80; 200; 600; 2000 nmol/l	6	pcs	blue (C1 - colourless)	2 months
3	CONTROL Control serum (0.8 ml)	1	pcs	colourless	2 months
4	CONJ HRP Conjugate, 11 ml	1	pcs	red	until exp. date
5	SUBS TMB Substrate solution, 11 ml	1	pcs	colourless	until exp. date
6	BUF WASH 21X Washing solution concentrate 21x, 22 ml	1	pcs	colourless	Concentrate - until exp. date Diluted washing solution - 1 month at 2...+8 °C or 5 days at RT
7	STOP Stop solution, 11 ml	1	pcs	colourless	until exp. date
8	N003 Plate sealing tape	2	pcs		N/A
9	K210I Instruction Cortisol EIA	1	pcs		N/A
10	K210Q QC data sheet Cortisol EIA	1	pcs		N/A

**5.2. Equipment and material required but not provided**

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 25–250 µl;
- Dry thermostat for +37 °C ±0,1 °C (optional: with shaking 600 rpm);
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0.

**5.3. Storage and stability of the Kit**

Store the whole kit at +2...+8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells **TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED)** to minimize exposure to moisture.

**6. SPECIMEN COLLECTION AND STORAGE**

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at +2...+8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

**7. TEST PROCEDURE****7.1. Reagent Preparation**

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 21X by 21 dilutions in distilled water.

**7.2. Procedural Note:**

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

**7.3. Assay flowchart**

See the example of calibration graphic in Quality Control data sheet.

**7.4. Assay procedure**

1	Put the desired number of microstrips into the frame; allocate 14 wells for the calibrators CAL 1-6 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Pipet 25 µl of calibrators CAL 1-6, control samples CONTROL and unknown samples into the wells.
3	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape (included into the kit).
4	Incubate 60 minutes at 37 °C .
5	Prepare washing solution by 21x dilution of washing solution concentrate (BUF WASH 21X) by distilled water. Wash the strips 5 times.
6	Dispense 100 µl of SUBS TMB into the wells
7	Incubate 10-20 minutes at +18...+25 °C
8	Dispense 100 µl of STOP into the wells.
9	Measure OD (optical density) at 450 nm.
10	Set photometer blank on air
11	Apply lin-log method for data reduction.

**7.5. Handling notes**

Calibrators and control sample(s) - only one freezing/thawing cycle is allowed  
 Incubation at +37 °C may be replaced by incubation at +37 °C with continuous shaking (ca 600 rpm) during 30 min

**8. QUALITY CONTROL**

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer’s instructions for use. Moreover the user must strictly adhere to the rules of *GLP* (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

**9. CALCULATION OF RESULTS**


**9.1.** Calculate the mean absorbance values (OD450) for each pair of calibrators and samples.

**9.2.** Plot a calibration curve on graph paper: OD versus cortisol concentration.

**9.3.** Determine the corresponding concentration of cortisol in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

**9.4.** Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

Calibrators	Value	Absorbance Units (450 nm) x 1000
CAL 1	0 nmol/l	3120
CAL 2	40 nmol/l	2092
CAL 3	80 nmol/l	1440
CAL 4	200 nmol/l	780
CAL 5	600 nmol/l	335
CAL 6	2000 nmol/l	138

**10. EXPECTED VALUES**

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for Cortisol. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units, nmol/l	
	Lower limit	Upper limit
Healthy donors	140	600

## 11. PERFORMANCE CHARACTERISTICS

### 11.1. Analytical specificity / Cross reactivity

Analyte	Cross-reactivity, % wt/wt
Cortisol	100
11-Deoxycortisol	0,9
Prednisolone	5,6
Corticosterone	0,6
11-Deoxycorticosterone	<0,1
Progesterone	<0,1
17-Hydroxyprogesterone	<0,1
Testosterone, Estradiol, Estriol	<0,1
Danazol	<0,01

### 11.2. Analytical sensitivity

Sensitivity of the assay was assessed as being 12 nmol/l.

### 11.3. Linearity

Linearity was checked by assaying dilution series of 5 samples with different cortisol concentrations. Linearity percentages obtained ranged within 90 to 110%.

### 11.4. Recovery

Recovery was estimated by assaying 5 mixed samples with known cortisol concentrations. The recovery percentages ranged from 90 to 110%.

## 12. LITERATURE

1. L. Thomas, Labor und Diagnose, 4. Auflage, 1992
2. Tietz, N.W., Textbook of Clinical Chemistry, Saunders, 1968
3. Foster, L. B. and Dunn, R.T. Clin. Chem: 20/3, 365 (1974)
4. De Lacerda, L., Kowarski, A., and Migeon, C.J. J. Clin. Endocr. and Metab: 36, 227 (1973)
5. Rolleri, E., Zannino, M., Orlandini, S., Malvano, R. Clin chim Acta 66 319 (1976)
6. Kobayashi, Y., et al Steroids, 32 no. 1 (1978)
7. Arakawa, H., Maeda, M., Tsuji, A. Anal. Biochem. 97 248 (1979)