Instruction for use

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**A SOLID-PHASE ENZYME IMMUNOASSAY  
FOR THE QUANTITATIVE DETERMINATION  
OF FREE PSA IN HUMAN BLOOD SERUM OR PLASMA**

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
**CONTENT**

1. INTENDED USE	9
2. SUMMARY AND EXPLANATION	9
3. PRINCIPLE OF THE TEST	10
4. WARNINGS AND PRECAUTIONS	10
5. KIT COMPONENTS	11
6. SPECIMEN COLLECTION AND STORAGE	12
7. TEST PROCEDURE	12
8. QUALITY CONTROL	14
9. CALCULATION OF RESULTS	14
10. EXPECTED VALUES	14
11. PERFORMANCE CHARACTERISTICS	15
12. LITERATURE	15

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1. Christensson A, Bjork T, Nilsson O, et al. Serum Prostate Specific Antigen Complexed to  $\alpha$ 1-Antichymotrypsin As An Indicator of Prostate Cancer. J. of Urol. 150:100-105; 1993.
2. Lilja H, Christensson A, Dahlen U, et al. Prostate-specific antigen in serum occurs predominantly in complex with  $\alpha$ -1-antichymotrypsin. Clin Chem 37:1618-1625, 1991.
3. Stenman U-H, Leinonen J, Alfthan H, Rannikko S, Tuhkanen K, Alfthan O. A complex between prostate-specific antigen and  $\alpha$ -1-antichymotrypsin is the major form of prostate-specific antigen in serum of patients with prostate cancer: assay of the complex improves clinical sensitivity for cancer. Cancer Res. 51:222-226, 1991.
4. Catalona WJ, Smith DS, Ratliff TL, Basler JW. Detection of organ-confined prostate cancer is increased through prostate-specific antigen-based screening. JAMA. 270:948-954, 1993.
5. Stamey TA, Yang N Hay AR, McNeal JE, Freiha, FS, Redwine E. Prostate-specific antigen as a serum marker for adenocarcinoma of the prostate. N Engl J Med. 317-: 909-916, 1987.
6. Junker R, Brandt B, Zechel C, and Assmann, G. Comparison of Prostate-specific antigen (PSA) measured by four combinations of free-PSA and total PSA assays. Clinical Chemistry 43:1588-1594, 1997.
7. Luderer AA, Chen Y-T, Soriano TF, et al. Measurements of the proportion of free to total prostate-Specific antigen improves diagnostic performance of prostate-specific antigen in the diagnostic gray zone of total prostate-specific antigen. Urology 46:187-194, 1995

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*Instruction for use*

## **A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF FREE PSA IN HUMAN BLOOD SERUM OR PLASMA**

### **1. INTENDED USE**

A solid-phase enzyme immunoassay for the quantitative determination of free PSA in blood serum or plasma.

This kit is designed for measurement of free PSA in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 42 unknown samples in duplicates.

### **2. SUMMARY AND EXPLANATION**

Prostate specific antigen (PSA) is a serin-like protease with molecular weight ca. 34 kDa and was initially found exclusively in normal prostatic gland as well as in prostatic fluid and seminal plasma. Later it was localized also in breast milk and, according to its enzymological properties, was classified as human prekallikrein 3.

In human serum, most of PSA forms complexes with serine protease inhibitor proteins (mostly alpha-1-antichymotrypsin, and alpha-2-macroglobulin). A minor proportion of PSA (free PSA) is circulating outside these complexes.

In a present test system, monoclonal antibodies PS2 capture both free and bound forms of PSA with an equal affinity ("equimolar" binding"). To detect captured free form of PSA, we use labelled monoclonal antibody PS1 which is highly specific for free PSA. The specificities and epitope mapping of these two antibodies were confirmed by independent laboratory (University of Turku, Finland).

It is generally accepted now that the ratio free PSA/ total PSA may help in differential diagnosis between benign diseases of prostatic gland (adenomas, inflammatory diseases) and carcinomas. This ratio should be measured when the total PSA level exceeds the used populational threshold level (usually, 4 ng/ml). In sera with lower concentrations of total PSA, especially in women, the correct measurement of free to total PSA ratio is not recommended by present test system.

### 3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to human PSA-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Unbound material is removed by washing procedure. Second antibodies – murine monoclonal to human free PSA, labelled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

### 4. WARNINGS AND PRECAUTIONS

**4.1.** For professional use only.

**4.2.** This kit is intended for in vitro diagnostic use only.

**4.3. INFECTION HAZARD:** There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

**4.4.** Avoid contact with stop solution containing 5,0%  $H_2SO_4$ . It may cause skin irritation and burns.

**4.5.** Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

**4.6.** Do not use the kit beyond the expiration date.

**4.7.** All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

**4.8.** Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

**4.9.** Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

**4.10.** Do not mix reagents from different lots.

**4.11.** Replace caps on reagents immediately. Do not swap caps.

**4.12.** Do not pipette reagents by mouth.

**4.13.** Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

**4.14.** Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

**4.15.** The Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

## 5.1. Contents of the Kit

## 5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1 SORB MTP	fPSA EIA strips, 8x12 wells	1	pcs		until exp.date
2 CAL 1-5	Calibrator set, 0.8 ml each. The set contains 5 calibrators: 0; 0.25; 0.75; 2.5; 5 ng/ml	5	pcs	blue (CI – colourless)	2 months
3 CONTROL	Control serum (0.8 ml)	1	pcs	colourless	2 months
4 CONJ HRP	Conjugate, 11 ml	1	pcs	purple	until exp.date
5 DIL	EIA buffer 11 ml	1	pcs	blue	until exp.date
6 SUBS TMB	Substrate solution, 11 ml	1	pcs	colourless	until exp.date
7 BUF WASH 21X	Washing solution concentrate 21X, 22 ml	1	pcs	colourless	Concentrate – until exp.date Diluted washing solution – 1 month at 2...+8 °C or 5 days at RT
8 STOP	Stop solution, 11 ml	1	pcs	colourless	until exp.date
9 N003	Plate sealing tape	2	pcs		N/A
10 K231I	Instruction fPSA EIA	1	pcs		N/A
11 K231Q	QC data sheet fPSA EIA	1	pcs		N/A

**5.2. Equipment and material required but not provided**

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 25–250 µl;
- Dry thermostat for +37 °C ±0.1 °C
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0.

**5.3. Storage and stability of the Kit**

Store the whole kit at +2...+8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells **TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED)** to minimize exposure to moisture.

**6. SPECIMEN COLLECTION AND STORAGE**

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at +2...+8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

**7. TEST PROCEDURE****7.1. Reagent Preparation**

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 21X by 21 dilutions in distilled water.

**7.2. Procedural Note:**

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

**7.3. Assay flowchart**

See the example of calibration graphic in Quality Control data sheet.

**7.4. Assay procedure**

1	Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1-5 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Pipet 50 µl of EIA buffer into each well.
3	Pipet 50 µl of calibrators CAL 1-5, control samples CONTROL and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
4	Incubate 30 minutes at 37 °C
5	Prepare washing solution by 21x dilution of washing solution concentrate BUF WASH 21X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times
6	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
7	Incubate 30 minutes at 37 °C
8	Wash the strips 5 times.
9	Dispense 100 µl of SUBS TMB into the wells
10	Incubate 10-20 minutes at +18...+25 °C
11	Dispense 100 µl of STOP into the wells.
12	Measure OD (optical density) at 450 nm.
13	Set photometer blank on first calibrator
14	Apply point-by-point method for data reduction.

**7.5. Handling notes**

Calibrators and control sample(s) - only one freezing/thawing cycle is allowed

## 8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

## 9. CALCULATION OF RESULTS


**9.1.** Calculate the mean absorbance values (OD450) for each pair of calibrators and samples.

**9.2.** Plot a calibration curve on graph paper: OD versus free PSA concentration.

**9.3.** Determine the corresponding concentration of free PSA in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

**9.4.** Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

Calibrators	Value	Absorbance Units (450 nm)
CAL 1	0 ng/ml	0.08
CAL 2	0.25 ng/ml	0.19
CAL 3	0.75 ng/ml	0.40
CAL 4	2.5 ng/ml	1.42
CAL 5	5 ng/ml	2.90

## 10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for fPSA. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Valid only for tPSA concentrations greater than 4.0 ng/ml

Ratio fPSA/tPSA (cat#K221, XEMA), %

Benign prostatic diseases	>10.0
Prostatic carcinoma	<10.0

## 11. PERFORMANCE CHARACTERISTICS

### 11.1. Analytical sensitivity

Sensitivity of the assay was assessed as being 0.07 ng/ml.

### 11.2. Linearity

Linearity was checked by assaying dilution series of 5 samples with different free PSA concentrations. Linearity percentages obtained ranged within 90 to 110%.

### 11.3. Recovery

Recovery was estimated by assaying 5 mixed samples with known free PSA concentrations. The recovery percentages ranged from 90 to 110%.

## 12. LITERATURE

1. Christensson A, Bjork T, Nilsson O, et al. Serum Prostate Specific Antigen Complexed to  $\alpha$ 1-Antichymotrypsin As An Indicator of Prostate Cancer. *J. of Urol.* 150:100-105; 1993.
2. Lilja H, Christensson A, Dahlen U, et al. Prostate-specific antigen in serum occurs predominantly in complex with  $\alpha$ -1-antichymotrypsin. *Clin Chem* 37:1618-1625, 1991.
3. Stenman U-H, Leinonen J, Alfthan H, Rannikko S, Tuhkanen K, Alfthan O. A complex between prostate-specific antigen and  $\alpha$ -1-antichymotrypsin is the major form of prostate-specific antigen in serum of patients with prostate cancer: assay of the complex improves clinical sensitivity for cancer. *Cancer Res.* 51:222-226, 1991.
4. Catalona WJ, Smith DS, Ratliff TL, Basler JW. Detection of organ-confined prostate cancer is increased through prostate-specific antigen-based screening. *JAMA.* 270:948-954, 1993.
5. Stamey TA, Yang N Hay AR, McNeal JE, Freiha, FS, Redwine E. Prostate-specific antigen as a serum marker for adenocarcinoma of the prostate. *N Engl J Med.* 317-: 909-916, 1987.
6. Junker R, Brandt B, Zechel C, and Assmann, G. Comparison of Prostate-specific antigen (PSA) measured by four combinations of free-PSA and total PSA assays. *Clinical Chemistry* 43:1588-1594, 1997.
7. Luderer AA, Chen Y-T, Soriano TF, et al. Measurements of the proportion of free to total prostate-specific antigen improves diagnostic performance of prostate-specific antigen in the diagnostic gray zone of total prostate-specific antigen. *Urology* 46:187-194, 1995

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Russian Diagnostic Manufacturers Association

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